


ФОТОГРАФИЧЕСКИЕ ИЗОБРАЖЕНИЯ МЕДИЦИНСКОГО ИЗДЕЛИЯ

**«Стент внутрисосудистый
Promus Element Plus»**

...на элемент 115 в описании...





Ультразвуковой Promus Element Plus в общем виде





систем, рентгеноконтрастные маркеры








Колпачок




Защитная трубка





Устройство защиты стента



Мандрен



PROMUS Element™ Plus

MONORAIL™

Everolimus-Eluting Coronary Stent System

atm - kPa Pressure		4.00mm Stent I.D. Stent O.D.	
11 - 1117	NOMINAL	3.96	4.11
12 - 1213		4.02	4.18
13 - 1317		4.08	4.23
14 - 1420		4.13	4.28
15 - 1517		4.17	4.32
16 - 1620	RATED*	4.21	4.36
17 - 1724		4.25	4.41
18 - 1827		4.30	4.46
19 - 1924		4.36	4.52
20 - 2027		4.43	4.58
21 - 2130		--	--
22 - 2227		--	--
*Rated Burst Pressure. DO NOT EXCEED.			

Internal BSC Bar Code



90375893-01



90634892-06D

**Boston
Scientific**


4.00mm
x8mm

PROMUS Element™ Plus
MONORAIL™

Everolimus-Eluting Coronary Stent System

[illegible]

atm. xPa	4.6mm	
Pressure	Stand D	Stand G
11 1113	NOMINAL 3.96	4.11
12 1213	4.02	4.18
13 1317	4.08	4.23
14 1420	4.13	4.28
15 1517	4.17	4.32
16 1620	RATED* 4.21	4.38
17 1724	4.25	4.41
18 1827	4.30	4.46
19 1924	4.36	4.52
20 2027	4.43	4.58
21 2130		
22 2227		



 $\leq 0.014\text{in (0.36mm)}$
 Recommended Guidewire
 $\geq 0.056\text{in (1.42mm)}$
 Recommended Guide Catheter
 (1) Includes Flushing Needle
 with Luer Fitting

 Contents (1)

GTIN

08714729795612

REF

Catalog No.

H7493918408400

STERILE/EO

Sterilized using
ethylene oxide

30 °C
15 °C
25 °C

Store at 25 °C (77 °F); excursions permitted to 15-30 °C (59-86 °F).

Do not open foil pouch until ready for use.

RH Protect from Humidity

PL Protect from Light

 See foil pouch and carton label for Use By date.

Internal BSC Bar Code

90634889-060

C € 0344

Made in IRELAND
Ballybrit Business Park
Galway, IRELAND



Boston
Scientific

PROMUS Element™ Plus

MONORAIL™

Everolimus-Eluting Coronary Stent System

Directions for Use	2
Instrucciones de uso	11
Mode d'emploi	21
Gebrauchsanweisung	31
Istruzioni per l'uso	41
Gebruiksaanwijzing	51
Instruções de Utilização	61

Boston
Scientific


t Cr

PROMUS Element™ Plus MONORAIL™

Everolimus-Eluting Coronary Stent System

Sistema de stent coronario con liberación de Everolimus, Système de stent coronaire à elution d'everolimus, Everolimus-eluting Coronary Stent System, Sistema di stent coronarico a rilascio di everolimus, Everolimus-eluierend coronaire stentsysteem, エベロリムス溶出型冠動脈ステントシステム, Everolimus-eluierende koronarstentsysteem, L'ustrojstvo elujiveho stentu excludující everolimus, Sistema de Stent Coronario para Elución de Everolimus, Everolimus-eluierende koronarstentsysteem, Everolimus-kibocsátó szirkizálójer-stentrendszer, System koronarnich stentu uvolňující everolimus, System stentu wierzchowego uwalniającego everolimus, Everolimus-eluierende koronarstentsysteem, Everolimus-Salinierte Koronar Stent System

Contents: (1)



GW ≤ 0.014 in (0.36mm)
Recommended Guidewire
GC ≥ 0.056 in (1.42mm)
Recommended Guide Catheter
NE (1) Includes Flushing Needle
with Luer Fitting

atm. kPa	Pressure	Stent I.D.	Stent O.D.
11	1117	NOMINAL 2.95	3.14
12	1213	3.01	3.20
13	1317	3.06	3.25
14	1420	3.10	3.29
15	1517	3.14	3.33
16	1620	RATED* 3.17	3.36
17	1724	3.20	3.39
18	1827	3.24	3.43
19	1924	3.28	3.47
20	2021	3.32	3.51
21	2130		
22	2221		

*Rated Burst Pressure. DO NOT EXCEED.

GTIN

08714729795599

LOT

19559689

REF

Catalog No.

H7493918408300

Use By

2018-08-02

30 °C
25 °C
15 °C

Store at 25 °C (77 °F);
excursions permitted to
15-30 °C (59-86 °F).

Do not open foil pouch
until ready for use

Protect from
Humidity

Protect from
Light

See foil pouch and cartons
label for Use By date



Картонная коробка – вид сбоку

