



ФОТО и Маркировка

«Контроль серии DxH 500 (DxH 500 Series Control)»


2017 г.

Общество с ограниченной ответственностью «Бекмен Культер»


109004, г. Москва, ул. Станиславского, д. 21, стр. 3
тел.: (495) 228-67-00, факс: (495) 228-67-01

ЭЛЕКТРОННАЯ КОПИЯ ОФОРМЛЕНА ЦЕНТРОМ ИСПЫТАНИЙ И РЕГИСТРАЦИИ МЕДИЦИНСКИХ ИЗДЕЛИЙ

E-mail: beckman.ru@beckman.com
www.beckmancoulter.com









[IVD]

Rx Only in the U.S.A.

Contents:

2 x 2.3 mL ABNORMAL LOW

2 x 2.3 mL NORMAL

2 x 2.3 mL ABNORMAL HIGH

2°C-8°C

DxH 500 Series Control is an assayed whole blood control designed to monitor values on multi parameter hematology cell counters. Refer to the Table of Expected Results for specific instrument models. / El control de la serie DxH 500 es un control de sangre entera analizada diseñado para controlar los valores en contadores celulares hematológicos de varios parámetros. Consulte los modelos de instrumentos específicos en la tabla de resultados esperados. / DxH 500 Serisi Kontrol, çok parametrelili hematoloji hücre saygıclarındaki değerleri izlemek için tasarlanan test edimiş tam kan kontrolüdür. Belirli cihaz modelleri için Beklenen Sonuçlar Tablosu'na bakın.

Package Insert/ Notice/ Gebrauchsinformation/ Foglieito illustrativo/ Prospecto/ Folhetos Informativos/ Pakkeindlæg/ Bopakstedel/ Pakkewedlegg/ Ένθετο ονομαστικό/ 添付文書/ 包装说明书/ Pakuotes įdėjimas/ Termidokijokortain/ Uloeka w opakowaniu/ Příložené informace/ Příbalový leták/ 첨부 설명서/ Ambalaj Prospekti/ Прокто-инструкци в упаковке/ Проктоеве/ 包裝說明書/ Complete pentru pachet/ Upulstvo u pakovanju/ Folheto informativo

Reactive Ingredients/ 成分:

Stabilized Human Blood/ 安定ヒト血液

Product contains biologically sourced material:
Human, Bovine and Porcine



WARNING

H315, H316, H319

Potassium Hydroxide < 1%

DxH 500 Series
Control



Manufactured for:
Beckman Coulter Ireland Inc.
Lismeehan
O'Callaghan's Mills
Co. Clare, Ireland



SDS



製造販売元: ベックマン・コールター株式会社
東京都江東区有明三丁目5番7号
TDC有明ウエストタワー

www.beckmancoulter.com/ru

LOT 160791 (01)11099190671014
2017-02-05 (17)170201
(10)160791



4H6280360721Q



4553170205160791QF




REF B36872




Printed in USA Made in USA
© 2016 Beckman Coulter, Inc.
B36871-AB



**BECKMAN
COULTER**



6



DxH 500 Series
Control



Manufactured for:
Beckman Coulter Ireland Inc.
Lismeehan
O'Callaghan's Mills
Co. Clare, Ireland

製造販売元: ベンタマシ・コールター株式会社
東京都江東区有明一丁目6番7号
TOKYO有明タニシトタワー

SDS



www.beckmancoulter.com/ifu

IVD

Rx Only In the U.S.A.

Contents:
2 x 2.3 mL, ABNORMAL LOW
2 x 2.3 mL, NORMAL
2 x 2.3 mL, ABNORMAL HIGH

DxH 500 Series Control is an assayed whole blood control designed to monitor values on multi parameter hematology cell counters. Refer to the Table of Expected Results for specific instrument models. / El control de la serie DxH 500 es un control de sangre entera analizada diseñado para controlar los valores en contadores celulares hematológicos de varios parámetros. Consulte los modelos de instrumentos específicos en la tabla de resultados esperados. / DxH 500 Serisi Kontrol, çok parametrelili hematoloji hücre saygınlardaki değerleri izlemek için tasarlanan test edimiş tam kan kontrolüdür. Belirli cihaz modelleri için Beklenen Sonuçlar Tablosu'na bakın.

Package Insert/ Notice/ Gebrauchsinformation/ Foglietto illustrativo/ Prospekt/ Folhetos informativos/ Pakkeindlæg/ Bipskedsedel/ Pakkevedlegg/ Evénro ovtovuatodig/ 添付文書/ 包装説明書/ Pakuotes įdėkės/ Termékdokumentáció/ Ulotka w opakowaniu/ Příložený listek/ 封入書/ 包装説明書/ Ambalaj Prospektüsü/ Тартов-сонзакшу с утасовоу/ Тартовса/ 包装内附説明/ Complete pentru pachet/ Uputstvo u pakovanju/ Folheto informativo

Reactive Ingredients/ 成分:
Stabilized Human Blood/ 安定ヒト血液
Product contains biologically sourced material:
Human, Bovine and Porcine



WARNING
H315, H316, H319
Potassium Hydroxide < 1%



REF B36872



Printed in USA Made in USA
© 2016 Beckman Coulter, Inc.
B36871-AB



HOLD AREA
3-7/8" x 1-3/8"

DxH 500 Series
Control

THIS BOX IS FOR TECHNICAL INFORMATIONAL USE ONLY. Refer to "4print.zip" file for all printing specifications.

Special Instructions: N/A

Description: CTN; DXH 500 SERIES CONTROL

Ink Color(s): 100% Black, PMS 285C, PMS 186C

Part Number: B36871

Revision: AB


Label Size: 3-7/8" x 2-3/8" x 4-1/2"

Scale: 1:1

IVD**2.3 mL ABNORMAL LOW****WARNING**
H315, H316, H319**LOT**2" x 7/16"
Hold Area for Barcode© 2016 Beckman Coulter, Inc.
Printed in USA B36875-AB**DxH 500 Series Control**

THIS BOX IS FOR TECH COMM INFORMATIONAL USE ONLY. Refer to "4print.zip" file for all printing specifications.

Special Instructions: Reduce 50%**Description:** LBL; PRI, DXH 500 SERIES CONTROL - LOW**Ink Color(s):** 100% Black, PMS 186C, PMS 2925C**Part Number:** B36875**Revision:** AB**Label Size:** 2" x 1-1/4"**Scale:** 2 : 1

IVD	WARNING H315, H316, H319				
2.3 mL NORMAL					
LOT					
2" x 7/16" Hold Area for Barcode					
© 2016 Beckman Coulter, Inc. Printed in USA B36877-AB					



DxH 500 Series Control

THIS BOX IS FOR TECH COMM INFORMATIONAL USE ONLY. Refer to "4print.zip" file for all printing specifications.			
Special Instructions: Reduce 50%			
Description: LBL; PRI, DXH 500 SERIES CONTROL - NORMAL		Ink Color(s): 100% Black, PMS 186C	
Part Number: B36877	Revision: AB	Label Size: 2" x 1-1/4"	Scale: 2 : 1

IVD

2.3 mL ABNORMAL HIGH

WARNING
H315, H316, H319



LOT



2" x 7/16"
Hold Area for Barcode

© 2016 Beckman Coulter, Inc.
Printed in USA B36874-AB



DxH 500 Series Control

THIS BOX IS FOR TECH COMM INFORMATIONAL USE ONLY. Refer to "4print.zip" file for all printing specifications.

Special Instructions: Reduce 50%

Description: LBL; PRI, DXH 500 SERIES CONTROL - HIGH

Ink Color(s): 100% Black, PMS 186C, PMS 356C

Part Number: B36874

Revision: AB

Label Size: 2" x 1-1/4"

Scale: 2 : 1

ПРОШито, пронумеровано
и скреплено печатью

13



Книжка листов
Иванова И.А.