

Пленка рентгеновская КХ, различных форматов, вариант исполнения КХ 123.

LUCKY KX123 MEDICAL X-RAY FILM

SPEC: (14 X 14) in. 35cm X 35cm QTY.:100 Sheets

- Storage**
The film should be kept in the original package and stored in dry, cool and no dust room without sunlight, harmful gases and radiation.
Details as follow:
1. The suitable temperature is 15°C-25°C, the humidity is 30%-65%/RH, and the film should be stored far from heat source.
2. Stored in room without hydro sulfuric gas, ammonia, sulfur dioxide and formaldehyde gas, as well as developing chemicals.
3. For heat-radiations and the radiation dose should be no more than 5.1 r-Dyph. It can bear normal radiation inspections at stations and airports, but not more than three times.
4. Stored carefully to avoid any harmful effect caused by pressure.
In above environments, the product can be stored for 2 years.
Caution
1. The best conditions of usage should be: temperature 18°C-24°C, humidity 50%/RH-65%/RH.
2. The film stored at low temperature should be conditioned at room temperature for more than 4 hours before opening the package in a dark room for using in order to prevent condensation on film surface and degradation of the film.
3. Handle carefully to avoid defect and static electricity discharge caused by pressure, creasing, buckling and friction, and not to touch film with wet hands.
4. Before use the film, check the type and the specification of the matching intensifying screen, keep the intensifying screen clean and without scratch.
5. When dispose the waste film and chemical liquid after used this product, it is obligatory to confirm to the requirements of local laws and regulations.
Processing
This film is applicable for the High-Speed processors with processing time of 80s, 120s or 180s, etc. The former is recommended with Lucky processing chemical G-20 or G-21, and other equivalent chemical, as recommended, developing at temperature 30°C-35°C.
Please see instruction for more details.



Lucky Healthcare Company Limited
Address: No. 8, Long South Street, Beijing, China
Tel: +86-10-53888881 Fax: +86-10-53888882
E-mail: lucky@lucky.com



Shanghai International Holding Corp GmbH(Europe)
Add: Effentasse 80, 21517 Hamburg, Germany
Tel: +49-40-2017170 Fax: +49-40-2017172
E-mail: shiholding@shiner.com



ISO 9001
ISO 13485

Пленка рентгеновская КХ, различных форматов, вариант исполнения КХ 171.

LUCKY KX171 MEDICAL X-RAY FILM

SPEC.: (12 X 16) in. 30cm X 40cm QTY.:100 Sheets

Storage

The film should be kept in the original package and stored in dry, cool and no dust rooms without sunlight, harmful gases and radiations.

Details as follow:

1. The storage temperature is 15°C-25°C, the humidity is 50%-65% RH, and the film should be stored for three years at most.
2. Except in cases where special instructions are given, the film should be stored in its original packaging in a dry, cool and no dust room.
3. Far from radiations and the radiation dose should be no more than 0.1 µSv/h. It can bear normal radiation inspection of shipment and storage, but not more than three times.
4. Stored vertically to avoid any harmful effect caused by pressure.

In above environments, the product can be stored for 3 years.

Caution

1. The best condition of usage should be temperature 15°C-24°C, humidity 50%-65% RH.
2. The film should be kept in its original packaging in a dry, cool and no dust room for more than 4 hours before opening the package in a dark room for using in order to prevent condensation on film surface and degradation of the film.
3. Handle carefully to avoid defect and static electricity discharge caused by pressure, creasing, buckling and friction, and not to touch film with wet hands.
4. Before use the film, check the type and the specification of the matching intensifying screen, keep the intensifying screen clean and without scratch.
5. When dispose the waste film and chemical liquid after used this product, it is obligatory to confirm to the requirements of local laws and regulations.

Processing

This film is applicable for the High-Speed processors with processing time of 90s, 120s or 150s, etc. The former is recommended with Lucky processing chemical G-30 or G-31, and other equivalent chemical as recommended, developing at temperature 33°C - 35°C.

Please see Instruction for more details



Lucky Healthcare Company Limited

Address: No. 4, Lele Street, Beijing, Hebei, China
Tel: +86-312-3328111 Fax: Code 071854

Web: www.luckyfilm.com



Shanghai International Holding Corp. GmbH(Europe)

Address: Eifelstrasse 60 70537 Mannheim Germany
Tel: +49-62-7011110 Fax: +49-62-7011115

E-Mail: shanghai@luckyfilm.com



10-100110-17

Простуновано и скреплено печатью

2 (два)

Листов

[Signature]
(подпись)

" " "



2018 г.