


TREATMENT CENTER

Ergonomy and Productivity


Easy access to the equipment's resources combined with the professional's and patient's appropriate positioning are fundamental conditions for the dental surgeon productivity and quality of life, preventing repetitive-strain and, at the same time, reducing attendance time.



1. Back to zero
2. 3 programmable work positions
3. Light (4 intensities)
4. Backrest and seat movements

INTUITIVE AND MULTI-FUNCTION FOOTSWITCH

3 work positions with synchronized movements that can be programmed by the operator. Including a program for switching the dental light on at your preferred light intensity.



EXCLUSIVE FLEXTYPE ARM (WITH A PNEUMATIC LOCK)

It provides easy positioning and locking of the delivery unit regardless of the work position used. It meets the requirements of all bio-types.

AMBIDEXTROUS ARM OF THE DELIVERY UNIT

It enables movements of up to 150° for the delivery unit, thus making it possible to use both sides of the chair with important ergonomics gain.

150°

PNEUMATIC DRIVE BLOCK

Made in noble metal, it prompts the tips to work automatically when they are taken from their holders and the foot control is pushed.

