

Фотографические изображения общего вида на медицинское изделие

**Калибраторы для обеспечения правильности количественного
определения витамина B12 иммунохемилюминесцентным методом в
сыворотке и плазме крови человека на иммунохимических анализаторах
Alinity i "B12 Калибраторы (Alinity i B12 Calibrators)"**

(размером не менее 18 x 24 сантиметра)

Alinity i B12 Calibrators

 **Abbott**

			CONC pg/mL	CONC pmol/L
CAL A	1 x 3.0 mL	0	0	0
CAL B	1 x 3.0 mL	110	81	81
CAL C	1 x 3.0 mL	250	184	184
CAL D	1 x 3.0 mL	500	369	369
CAL E	1 x 3.0 mL	1000	738	738
CAL F	1 x 3.0 mL	2000	1476	1476



CONTAINS: AZIDE



www.abbottdiagnostics.com/IFU

REF

07P6701



Exp.

2020-03-15


LOT

03604UI00

R01




G70994R02



(01) 00380740131388 (17) 200315
(10) 03604UI00 (240) 07P6701 (91) B12





Alinity i B12 Calibrators

Created December 2017.

NAME

Alinity i B12 Calibrators (also referred to as B12 Cals)

INTENDED USE

The Alinity i B12 Calibrators are used to calibrate the Alinity i analyzer when used for the quantitative determination of vitamin B12 in human serum and plasma using the Alinity i B12 Reagent Kit. For additional information, refer to the Alinity i B12 reagent package insert or the Alinity i-Series Operations Manual.

The following warnings and precautions apply to: CAL A . CAL P	
DANGER	Contains disodium tetraborate decahydrate and sodium azide.
H360	May damage fertility or the unborn child.
EUH032	Contact with acids liberates very toxic gas.
Prevention	
P201	Obtain special instructions before use.
P273	Prevent release into the environment.
P501	Dispose of contents and container in accordance with local, state, and federal regulations.



Alinity i B12 Calibrators

		CONC pg/mL	CONC pmol/L
CAL A	1 x 3.0 mL	0	0
CAL B	1 x 3.0 mL	110	81
CAL C	1 x 3.0 mL	250	184
CAL D	1 x 3.0 mL	500	369
CAL E	1 x 3.0 mL	1000	738
CAL F	1 x 3.0 mL	2000	1476

CONTAINS: AZIDE

REF 07P6701

Exp. 2020-03-15

LOT 03604UI00

R01

G70994R02

2°C — 8°C

(01) 00380740131388 (17) 200315

(10) 03604UI00 (240) 07P6701 (91) B12



всего прошито, при
и скреплено печать
6
лист

