

Макет этикетки потребительской упаковки

	Набор реагентов для иммунохроматографического выявления антигена SARS-CoV-2 в мазке из носоглотки (NowCheck COVID-19 Ag Test). Lot № 1901D006. БиоНоут, Инк. (BioNote, Inc.) 22, Самсунг 1-ро 4-гил, Хвасон-си, Генджи-до, 18449 Республика Корея Телефон: +82-31-211-0516 Факс: +82-31-8003-0618 Интернет: www.bionote.co.kr Представитель производителя в РФ: ООО «Авивир» 141401, Московская обл., г. Химки, ул. Рабочая, д. 2А, стр. 1, этаж/пом. 2/7 Тел (495) 925-30-74	
Для профессионального применения Состав набора: 1. Тестовая кассета (в индивидуальной упаковке из фольги в комплекте с осушителем) – 25 шт.; 2. Пробирка с буфером для экстрагирования биологического образца – 25 шт.; 3. Насадка с капельницей – 25 шт.; 4. Стерильный тампон для взятия мазка – 25 шт.; 5. Подставка из картона – 1 шт.; 6. Набор пленочных этикеток – 1 шт.; 7. Инструкция к медицинскому изделию – 1 шт.		
REF	LOT 	
 PY	Дата регистрации 	
CE	     25   	



 NowCheck
COVID-19 Ag Test

 REF RG1901

 Lot No.: 1901D006, Code: GEN
 Mfg. Date: Sep.03.2020
 Exp. Date: Sep.02.2022

 
Autonizna Republika Slovenija
MT Promed Consulting GmbH
Avtorizirani za delo na terenu
Tel: +386 (0)4 681 11 111 Fax: +386 (0)4 681 11 112
E-mail: info@mt-promed.si

Test procedure

1. Insert a sterile swab into the center of the nostril, aiming for the septum.
2. Rotate the swab head to collect mucus.
3. Insert the swab into the extraction buffer tube, ensuring the tip is submerged.
4. Squeeze the tube to mix the sample.
5. Insert the swab into the test cassette.
6. Wait for 15 minutes.

Test result

Negative

Positive

Invalid

For more information, please visit www.nowcheck.com.
 Do not use the test cassette after the expiration date.
 Do not use the test cassette if the result is unclear.



EXTRACTION BUFFER TUBE

NowCheck COVID-19 Ag Test

IVD



1901D006



Sep.03,2020




Sep.02,2022

CE



BIONOTE

Regulatory Affairs





NOZZLE CAP

NowCheck COVID-19 Ag Test

LOT 1901D006
Sep.03,2020
Sep.02,2022

IVD
CE

BIONOTE (Republic of Korea)




Peel Here
Sterile only if peel is intact



iClean®

Specimen Collection Swab

STERILE EO



Product Code: 12-01
MFG: 20200510
EXP: 20270509
LOT: 202005001


Peel Here
Sterile only if peel is intact



iClean®

Specimen Collection Swab

STERILE EO



CE FDA

Product Code: 12-01
MFG: 20200510
EXP: 20270509
LOT: 202005001

Peel Here
Sterile only if peel is intact



iClean®


Specimen Collection Swab

STERILE EO




CE FDA

Product Code: 12-01
MFG: 20200510
EXP: 20270509
LOT: 202005001



Please attach this film on the result window of device, when performing the test outdoor or in an environment with lot of air circulation (clean bench).

L24RT2ENR0 / Issue date: 2020.04



Прошито, пронумеровано и
скреплено печатью 12 листов.

Генеральный директор
ООО «АВИВИР»

Рабинович Н.М.

