

Фотографические изображения медицинского изделия

«Наборы инфузионные Superline, EasyRelease, Soft-Release-ST и резервуар для инсулина DANA 3cc Syringe Insulin Reservoir, одноразовые стерильные, для помпы инсулиновой DANA Diabecare, в вариантах исполнения»

производства «СООИЛ Девелопмент Ко., Лтд.», Корея / «SOOIL Development Co., Ltd.», Korea




Фото 1. Набор инфузионный Superline
в индивидуальной блистерной упаковке и без упаковки.

SUPERLINE

For use only with
the DANA DIABECARE insulin pump.

Cat. No : SU 201

Rectangular Needle

CE
0197

Tubing length : 550mm(21.65inch)
Needle length : 6.5mm(0.25inch)
Inner diameter : 0.4mm(0.02inch)
Outer diameter : 1.5mm(0.06inch)
NHRIC Code : 8416.0201.30

CAUTION!

This device is restricted to sale by or on the order of a physician.
Change the reservoir and the infusion set regularly; at least every 72 hours.
Check blood glucose 1~2 hours after insertion.
If blood glucose is high, change the set.
Keep dry at room temperature.

STERILE EO

LOT

120223



Do not reuse



2016.02.22



2013.02.23

SOOIL

www.sooil.com

Made in Korea(South)

Manufactured by Sooil Development Co., Ltd.

111-1 Heukseok-dong, Dong jak-gu, Seoul, Korea / Tel: +82-2-3463-0041, Fax: +82-2-3463-7707

Distributed by

SOOIL, LLC: 5677 Oberlin Drive, Suite 101, San Diego, CA92121 / Tel: 1-866-747-6645

MDCI: Arundel House, 1 Liverpool Gardens, Worthing, West Sussex BN11 1SL, UK / Tel: 44(0)-293-429608

Фото 2. Набор инфузионный Superline. Маркировка индивидуальной блистерной упаковки.



Фото 3. Набор инфузионный Superline. Маркировка потребительской картонной упаковки.



Фото 4. Набор инфузионный Superline. Маркировка потребительской картонной упаковки.




Фото 5. Набор инфузионный EasyRelease
в индивидуальной блистерной упаковке и без упаковки.

EasyRelease

For use only with
the DANA DIABECARE insulin pump.

REF SU 313

Rectangular Needle

CE
0120

Tubing length : 700mm(27.55inch)
Needle length : 7mm(0.28inch)
Inner diameter : 0.4mm(0.02inch)
Outer diameter : 1.5mm(0.06inch)
NHRIC Code : 8416.313.30

CAUTION!

This device is restricted to sale by or on the order of a physician.
Change the reservoir and the infusion set regularly; at least every 72 hours.
Check blood glucose 1~2 hours after insertion.
If blood glucose is high, change the set.
Keep dry at room temperature.

STERILE EO

LOT

131025

  Do not reuse



2016. 10. 24



2013. 10. 25

SOOIL

www.sooil.com

Made in Korea(South)

Sooil Development Co., Ltd.

111-1, Heukseok-dong, Dongjak-gu, Seoul, Korea / Tel: 82(2)-3463-0041 / Fax: 82(2)-3463-7707

Distributed by

SOOIL, LLC: 5677 Oberlin Drive, Suite 101, San Diego, CA92121 / Tel: 1-866-747-6645

 RSQR Ltd., Ludgate House, 107-111 Fleet Street, London EC4A2AB

www.rsqa.co.uk, info@rsqa.co.uk, Tel: +44(0)207-936-9412

LA-UG-203.rev.06

Фото 6. Набор инфузионный EasyRelease. Маркировка индивидуальной блистерной упаковки.




Фото 7. Набор инфузионный EasyRelease. Маркировка потребительской картонной упаковки.




Фото 8. Набор инфузионный EasyRelease II
в индивидуальной блистерной упаковке и без упаковки.

EasyRelease II

For use only with
the DANA DIABECARE insulin pump.

Cat. No : SU 401

Rectangular Needle

CE
0120

Tubing length : 700mm(27inch)
Needle length : 7mm(0.28inch)
Inner diameter : 0.4mm(0.02inch)
Outer diameter : 1.5mm(0.06inch)
NHRIC Code : 8416.0401.30


CAUTION!

This device is restricted to sale by or on the order of a physician.
Change the reservoir and the infusion set regularly; at least every 72 hours.
Check blood glucose 1~2 hours after insertion.
If blood glucose is high, change the set.
Keep dry at room temperature.

STERILE EO

LOT

130823



Do not reuse



2016. 08. 22



2013. 08. 23

SOOIL

www.sooil.com

Made in Korea(South)

Sooil Development Co., Ltd.

111-1, Heukseok-dong, Dongjak-gu, Seoul, Korea / Tel: 82(2)-3463-0041 / Fax: 82(2)-3463-7707

Distributed by

SOOIL, LLC; 5677 Oberlin Drive, Suite 101, San Diego, CA 92121 / Tel: 1-866-747-6645

RSQR Ltd., Ludgate House, 107-111 Fleet Street, London EC4A 2AB

www.rsqa.co.uk, info@rsqa.co.uk, Tel: +44(0)207-936-9412

LA-UC-301 rev.05

Фото 9. Набор инфузионный EasyRelease II. Маркировка индивидуальной блистерной упаковки.




Фото 10. Набор инфузионный EasyRelease II. Маркировка потребительской картонной упаковки.




Фото 11. Набор инфузионный Soft-Release-ST
в индивидуальной блистерной упаковке и без упаковки

Soft-Release-ST

For use only with
the DANA DIABECARE insulin pump.

Cat. No : SR 102

Straight Needle

CE
0120

Tubing length : 700mm(27inch)
Needle length : 19mm(0.75inch)
Inner diameter : 0.4mm(0.02inch)
Outer diameter : 1.5mm(0.06inch)
NHRIC Code : 8416.0102.31

CAUTION!

This device is restricted to sale by or on the order of a physician.
Change the reservoir and the infusion set regularly; at least every 72 hours.
Check blood glucose 1~2 hours after insertion.
If blood glucose is high, change the set.
Keep dry at room temperature.

STERILE EO



Do not reuse

SOOIL

www.sooil.com

Made in Korea(South)

LOT

131212



2016. 12. 11



2013. 12. 12

Sooil Development Co., Ltd.

111-1, Heukseok-dong, Dongjak-gu, Seoul, Korea / Tel: 82(2)-3463-0041 / Fax: 82(2)-3463-7707

Distributed by

SOOIL, LLC: 5677 Oberlin Drive, Suite 101, San Diego, CA92121 / Tel: 1-866-747-6645

RSQR Ltd., Ludgate House, 107-111 Fleet Street, London EC4A2AB

www.rsqa.co.uk, info@rsqa.co.uk, Tel: +44(0)207-936-9412

LA-UG-401 rev. 03

Фото 12. Набор инфузионный Soft-Release-ST. Маркировка индивидуальной блистерной упаковки.

10 pcs

INFUSION SET FOR DANA DIABECARE INSULIN PUMP

Soft-Release-ST

10 pcs

CE
0120

Soft-Release-ST

INFUSION SET FOR DANA DIABECARE INSULIN PUMP

SINGLE USE ONLY, STERILIZED BY E.O. GAS

Manufactured by SOOIL Development Co., Ltd.
111-1, Heukseok-dong, Dongjak-gu, Seoul, Korea
Tel: 82(2)-3463-0041 / Fax: 82(2)-3463-7707

Distributed by
SOOIL, LLC: 5677 Oberlin Drive, Suite 101, San Diego, CA92121
Tel: 1-866-747-6645
1000 MOORE Avenue House 11 Cornwell Gardens, Worthing, West Sussex BN11 1SL, UK
Tel: 44(0) 1903 774777

CAUTION!

This Device is restricted to sale by or on the order of a physician.
Change the infusion set and catheter every 72 hours.
Check blood glucose 1-2 hours after insertion.
If blood glucose is high, change the set.
Keep dry at room temperature.

SOOIL Development Co., Ltd.

PA-UG-4
10 pcs

Soft-Release-ST

INFUSION SET FOR DANA DIABECARE INSULIN PUMP

SINGLE USE ONLY, STERILIZED BY E.O. GAS

Non-Toxic & Non-Pyrogenic

LOT
131212
2016.12.11
2013.12.12

SR102



Фото 13. Набор инфузионный Soft-Release-ST. Маркировка потребительской картонной упаковки.




Фото 14. Резервуар для инсулина DANA 3cc Syringe Insulin Reservoir
в индивидуальной блистерной упаковке и без упаковки

Reservoir Cap Linking Screw

DANA 3cc Syringe

22G X $\frac{1}{2}$ "

Caution!

- ▶ Do Not Use if Packaging Has Been Opened or Damaged
- ▶ For Use Only With Dana Diabecare Insulin Pump
- ▶ Change your syringe at least every 72 hours or as recommended by your health care professional.
- ▶ Keep Dry - Store at Room Temperature

Manufactured by

111-1, Heuksok-dong, Dongjak-gu, Seoul, Korea / Tel: 82(2)-3463-0041 / Fax: 82(2)-3463-7707

made in Korea

Distributed by

SDOIL, LLC: 5677 Oberlin Drive, Suite 101, San Diego, CA92121 / Tel: 1-866-747-6645

EC REP RSQR Ltd., Ludgate House, 107-111 Fleet Street, London EC4A2AB

www.rsqa.co.uk, info@rsqa.co.uk / Tel: +44(0)207-836-9412

Non-Toxicity
Non-Pyrogenic
Use U-100 Insulin only

STERILE EO

CE
0120



Single
use
only

LA-UG-701, Rev7

PEEL HERE

131024
2015.12.23
2013.12.24

Cal NO : DAC-002-01

LOT

PEEL HERE

Фото 15. Резервуар для инсулина DANA 3cc Syringe Insulin Reservoir.
Маркировка индивидуальной блистерной упаковки.

DANA 3cc Syringe



10 pcs

Reservoir : 3ml
Gauge : 22G X 1/2
Non-Toxicity &
Non-Pyrogenic

DANA 3cc Syringe

INSULIN RESERVOIR FOR DANA DIABECARE

- ▶ Do Not Use if Packaging Has Been Opened or Damaged
- ▶ For Use Only With Dana Diabecare Insulin Pump
- ▶ Change your syringe at least every 72 hours or as recommended by your health care professional.
- ▶ Keep Dry - Store at Room Temperature

LOT

3cc Syringe

Made in Korea (South)

121224

2016.12.23

2013.12.24



8 809220 587016

SOOIL Development Co.,Ltd

STERILE EO



SOOIL Development Co., Ltd.

제품명 : 주사기 품목허가번호 : 제 06-732호

제조원 : (주)수일개발

서울시 동작구 흑석동 111-1

TEL : (02)3463-0041, FAX : (02)3463-7707

Фото 16. Резервуар для инсулина DANA 3cc Syringe Insulin Reservoir.
Маркировка потребительской картонной упаковки.




Фото 17. Инфузионная система в сборе (набор инфузионный EasyRelease, помпа DANA Diabecare IIS, установленный резервуар для инсулина DANA Syringe Insulin Reservoir).

Прошито, пронумеровано
и скреплено печатью

18 лист (ов)

Генеральный директор
ООО «МедТехКонсалт»
И В. Копцева

