



МАРКИРОВКА МЕДИЦИНСКОГО ИЗДЕЛИЯ

Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1
с принадлежностями модели TEST1 THL/BCL/SDL/YDL/MDL

ВНУТРЕННЯЯ МАРКИРОВКА АНАЛИЗАТОРА

Модель, местоположение	Образец маркировки (оригинал)	Образец маркировки (русский язык)
а) Маркировка модели TEST1 THL на задней панели анализатора	<p>TEST1 THL REF SI 195.210/THL SN T1XXXX YYYY-MM-DD 115/230 V~ 50/60 Hz 225 VA (x2) T5A L 250V 5x20mm (for USA and Canada 115 VAC only) (pour USA et Canada seulement a 115CAV) 42 kg 49x54x60 cm Alifax S.r.l. via Merano, 30 - 33045 Nimis (UD) - ITALY CE IVD R ONLY (USA) Made in ITALY</p>	<p>Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 Каталожный номер/модель TEST1 SI 195.210/THL Серийный номер: T1XXXX ПУ №ФСЗ 2008/02628 от XXXX Электропитание 115/230 В, 50/60 Гц 225ВА (x2) T5A L 250V 5x20mm ДхШхВ 49x54x60 см Место производства: АЛИФАКС С.р.л. Вес нетто: 42 кг Алифакс С.р.л. Виа Петрарка, 2/1, 35020 Виа Мерано, 30, 33045 Польверара, Падуя, Италия Нимис (Удине), Италия Сделано в Италии Уполномоченный представитель изготовителя в РФ - ООО "Алифакс" 125367, г. Москва, ул. Габричевского, д.5, к. 1, комната 4-5, 21-22, 3-эт, пом. 1 +7(495)-544-50-55, info.russia@alifax.com ГТТГ-ММ-ДД</p>
б) Маркировка модели TEST1 BCL на задней панели анализатора	<p>TEST1 BCL REF SI 195.220/BCL SN T1XXXX YYYY-MM-DD 115/230 V~ 50/60 Hz 225 VA (x2) T5A L 250V 5x20mm (for USA and Canada 115 VAC only) (pour USA et Canada seulement a 115CAV) 47 kg 49x54x60 cm Alifax S.r.l. via Merano, 30 - 33045 Nimis (UD) - ITALY CE IVD R ONLY (USA) Made in ITALY</p>	<p>Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 Каталожный номер/модель TEST1 SI 195.220/BCL Серийный номер: T1XXXX ПУ №ФСЗ 2008/02628 от XXXX Электропитание 115/230 В, 50/60 Гц 225ВА (x2) T5A L 250V 5x20mm ДхШхВ 49x54x60 см Место производства: АЛИФАКС С.р.л. Вес нетто: 47 кг Алифакс С.р.л. Виа Петрарка, 2/1, 35020 Виа Мерано, 30, 33045 Польверара, Падуя, Италия Нимис (Удине), Италия Сделано в Италии Уполномоченный представитель изготовителя в РФ - ООО "Алифакс" 125367, г. Москва, ул. Габричевского, д.5, к. 1, комната 4-5, 21-22, 3-эт, пом. 1 +7(495)-544-50-55, info.russia@alifax.com ГТТГ-ММ-ДД</p>

с) Маркировка модели TEST1 SDL на задней панели анализатора

ALI FAX **TEST1 SDL** (01) 0 8056040 14152 6

REF SI 195.230/SDL

SN T1XXXX

YYYY-MM-DD

115/230 V~ 50/60 Hz 225 VA
(x2) T5A L 250V 5x20mm
(for USA and Canada 115 VAC only)
(pour USA et Canada seulement à 115CAV)

47 kg
Dimensions (LxWxH)
49x54x60 cm

AliFax S.r.l.
via Merano, 30 - 33045
Nimis (UD) - ITALY

CE **IVD** **Rx ONLY (USA)**
Made in ITALY

ALI FAX Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1

REF Каталогный номер/модель TEST1 SI 195.230/SDL

SN Серийный номер: T1XXXX

ПУ №ФСЗ 2008/02628 от XXXX
Электропитание 115/230 В, 50/60 Гц 225ВА
(x2) T5A L 250V 5x20mm ДхШхВ 49x54x60 см

Место производства: АЛИФАКС С.р.л. Вес нетто: 47 кг
Алифакс С.р.л. Виа Петрарка, 2/1, 35020
Виа Мерано, 30, 33045 Польверара, Падуа, Италия
Нимис (Удине), Италия Сделано в Италии

Уполномоченный представитель изготовителя в РФ - ООО "Алифакс"
125367, г. Москва, ул. Габричевского, д.5, к. 1, комната 4-5, 21-22, 3-эт.
пом. 1 +7(495)-544-50-55 info.russia@alifax.com

CE **IVD** ГТТТ-ММ-ДД

д) Маркировка модели TEST1 YDL на задней панели анализатора

ALI FAX **TEST1 YDL** (01) 0 8056040 14032 1

REF SI 195.240/YDL

SN T1XXXX

YYYY-MM-DD

115/230 V~ 50/60 Hz 225 VA
(x2) T5A L 250V 5x20mm
(for USA and Canada 115 VAC only)
(pour USA et Canada seulement à 115CAV)

47 kg
Dimensions (LxWxH)
49x54x60 cm

AliFax S.r.l.
via Merano, 30 - 33045
Nimis (UD) - ITALY

CE **IVD** **Rx ONLY (USA)**
Made in ITALY

ALI FAX Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1

REF Каталогный номер/модель TEST1 SI 195.240/YDL

SN Серийный номер: T1XXXX

ПУ №ФСЗ 2008/02628 от XXXX
Электропитание 115/230 В, 50/60 Гц 225ВА
(x2) T5A L 250V 5x20mm ДхШхВ 49x54x60 см

Место производства: АЛИФАКС С.р.л. Вес нетто: 47 кг
Алифакс С.р.л. Виа Петрарка, 2/1, 35020
Виа Мерано, 30, 33045 Польверара, Падуа, Италия
Нимис (Удине), Италия Сделано в Италии

Уполномоченный представитель изготовителя в РФ - ООО "Алифакс"
125367, г. Москва, ул. Габричевского, д.5, к. 1, комната 4-5, 21-22, 3-эт.
пом. 1 +7(495)-544-50-55 info.russia@alifax.com

CE **IVD** ГТТТ-ММ-ДД

е) Маркировка модели TEST 1 MDL на задней панели анализатора

ALI FAX **TEST1 MDL** (01) 0 8056040 14269 1

REF SI 195.250/MDL

SN T1XXXX

YYYY-MM-DD

115/230 V~ 50/60 Hz 225 VA
(x2) T5A L 250V 5x20mm
(for USA and Canada 115 VAC only)
(pour USA et Canada seulement à 115CAV)

47 kg
Dimensions (LxWxH)
49x54x60 cm

AliFax S.r.l.
via Merano, 30 - 33045
Nimis (UD) - ITALY

CE **IVD** **Rx ONLY (USA)**
Made in ITALY

ALI FAX Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1

REF Каталогный номер/модель TEST1 SI 195.250/MDL

SN Серийный номер: T1XXXX

ПУ №ФСЗ 2008/02628 от XXXX
Электропитание 115/230 В, 50/60 Гц 225ВА
(x2) T5A L 250V 5x20mm ДхШхВ 49x54x60 см


Место производства: АЛИФАКС С.р.л. Вес нетто: 47 кг
Алифакс С.р.л. Виа Петрарка, 2/1, 35020
Виа Мерано, 30, 33045 Польверара, Падуа, Италия
Нимис (Удине), Италия Сделано в Италии

Уполномоченный представитель изготовителя в РФ - ООО "Алифакс"
125367, г. Москва, ул. Габричевского, д.5, к. 1, комната 4-5, 21-22, 3-эт.
пом. 1 +7(495)-544-50-55 info.russia@alifax.com

CE **IVD** ГТТТ-ММ-ДД

Передняя этикетка TEST1 THL на лицевой стороне анализатора, на пластиковой крышке	  Test 1 THL	
Передняя этикетка TEST1 BCL на лицевой стороне анализатора, на пластиковой крышке	  Test 1 BCL	
Передняя этикетка TEST1 SDL на лицевой стороне анализатора, на пластиковой крышке	  Test 1 SDL	
Передняя этикетка TEST1 YDL на лицевой стороне анализатора, на пластиковой крышке	  Test 1 YDL	
Передняя этикетка TEST1 MDL на лицевой стороне анализатора, на пластиковой крышке	  Test 1 MDL	

ВНЕШНЯЯ (ТРАНСПОРТНАЯ) МАРКИРОВКА (на примере модели TEST1 THL)

		  	
TEST1 THL		 	
REF	SI 195.210/THL	Only for US market	
SN	T1xxx	 65 °C / 149 °F -20 °C / -4 °F	
MODEL: [X]THL []BCL []SDL []SDL-S []YDL []MDL []XDL		 80% 0%	
FIRMWARE VERSION: _____		 	
DESTINATION: CUSTOMER		(01) 0 8056040 14151 9 (11) (21) T1xxx	
	 Alifax S.r.l. via Merano, 30 - 33045 NIMIS (UD) - ITALY	 YYYY-MM-DD	

ESR_EL_SI195210THL_TEST1THL_2.0

		  	
Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1, МОДЕЛЬ TEST1 THL			
REF	№ по кат. SI 195.210/THL	 65 °C -20 °C	
SN	Серийный номер T1XXX	 80 0	
РУ №ФСЗ 2008/02628 от XXXX года		 	
ПРЕДНАЗНАЧЕН ДЛЯ:   Алифакс С.р.л. Виа Мерано, 30-33045 НИМИС (УД) - ИТАЛИЯ		 ГГГГ-ММ-ДД	(01) 0 8056040 14151 9 (11) (21) T1XXX
Уполномоченный представитель производителя в РФ: ООО «АЛИФАКС», адрес: 125367, г. Москва, ул. Габричевского, д.5. к.1, комн.4-5, 21-22, 3-эт, пом. 1, телефон +7 (495) 544-50-5, info.russia@alifax.com			



ООО «АЛФАКС»

В данном документе
пронумеровано и прошито

4 (Четыре)
лист(а)ов



ООО «АЛФАКС»

**Фотографические изображения
медицинского изделия**

**Анализатор автоматический для
определения скорости оседания
эритроцитов (СОЭ) TEST1 с
принадлежностями**

(Регистрационное удостоверение

№ ФСЗ 2008/02628 от 17 июня 2021 года)



Фото 1 Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 THL (общий вид)



Фото 2 Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 THL (фото маркировки)



TEST 1

Руководство по эксплуатации

Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 с принадлежностями

Действительно для всех моделей

Применимо к версии 9.00X



Редакция для Российской Федерации



**Прибор для профессионального использования в области
In Vitro диагностики**

Copyright © Alifax S.r.l.

Информация в данном руководстве является собственностью ALIFAX. Все права защищены.
Несанкционированное копирование руководства или его частей запрещено.
Торговая марка TEST 1 является собственностью Alifax S.r.l.



Фото 4 Руководство по эксплуатации на электронном носителе (CD-диск)



Фото 5 Руководство по эксплуатации на электронном носителе (флеш-карта)



FASCICOLO APPARECCHIO D195.08.03.02 Rev 1.9

DEVICE HISTORY RECORD

Pag. 1/1

TEST1 THL / BCL / SDL / YDL / MDL / XDL / HDL (SI 195.2XX)

s.n. 713232commessa (Job Order) 2023-211280Data PRD 22/2/23

CORREDO DI SERIE / PACKING LIST

STRUMENTO INSTRUMENT	MODELLO MODEL						LINGUA LANGUAGE	
	THL	BCL	SDL	YDL	MDL	XDL	EN	IT
☒ TEST1 SI 195.2XX	☒	☐	☐	☐	☐	☐	☒	☐

Corredo di serie / Packing list

- ☒ N.° 2 SI195800 Rotoli Carta termica / Thermal Paper rolls
- ☒ N.° 2 SI205801 Tanica di lavaggio / Wash Tank
- ☒ N.° 2 SI205801 Tanica di scarico / Waste Tank
- ☒ N.° 2 Fusibili di alimentazione TSA / Main fuses TSA
- ☒ N.° 1 Cavo Alimentazione / Main cable plug
- ☐ Cavo alimentazione personalizzato/Personalized main cable _____
- ☒ N.° 1 Ago a vite / Needle
- ☒ N.° 1 CD Manuale operativo multilingua / CD Multi-Language Operative manual (Language and version)
ESR_UM_SI195_TEST1_SW 11-COA-4-0-EN-ES-IT-PT-FR-DE-FI ESR-CD-SI195-TEST1-SW-M-004-1-1
- ☒ ESR_IQOQ_SI195_TEST1: Modulo di installazione / Installation Report

SOLO STRUMENTI SDL / SDL INSTRUMENTS ONLY:

- ☐ N.° 4 Pettine rack Sysmex XN / Sysmex XN rack scallop
- ☐ Nota tecnica n.° 2004 EN / Technical note n.° 2004 EN

SOLO STRUMENTI SDL PER PROVETTE SARSTEDT / SDL INSTRUMENTS FOR SARSTEDT TUBES ONLY:

- ☐ N.° 4 Pettine x Sarstedt 11X66 / Sarstedt 11x66 scallop
- ☐ Nota tecnica n.° 2006 EN / Technical note n.° 2006 EN

- ☒ Dichiarazione di conformità CE / EC declaration of conformity
- ☒ DPRD0105: Prove elettriche di routine / Electrical routine tests
- ☒ Pag 8/9: Certificato di Collaudo / Final Test report
- ☒ Pag 9/9: Packing List / This document

Compilatore/Signature CassellaData/Date 16/03/2023

Approvato dal Responsabile Produzione Strumenti per rilascio finale / Manufacturer Manager Approval

DATA/Date 16/03/2023 FIRMA/Signature [Signature]

Approvato da Assicurazione Qualità (Approved by Quality Assurance)

Data 16.03.2023
(Date)Firma [Signature]
(Signature)

ALIFAX S.r.l.
Via Petrarca 2/1 - 35020 POLVERARA (PD) - ITALY
Tel. +39 0490992000 e-mail: info@alifax.com
PI 04337640280



ALIFAX S.r.l.
Via Merano 30 - 33045 NIMIS (UD) - ITALY
Documents available on: www.alifax.com

©Copyright 2022. Documento segreto non divulgabile a terzi, di proprietà esclusiva di ALIFAX S.r.l. Ogni diritto riservato.
©Copyright 2022. This document is exclusive property of ALIFAX S.r.l, confidential and not freely disclosable. All rights reserved.

Фото 6 Упаковочный лист (анализатор TEST1 модель THL)



INSTALLATION REPORT for TEST1 Instruments

Rev. 3.0 - Valid since 2022 May 25

Dear Distributor,

You are kindly requested to fill-in this form during the installation and return a copy within one month since the installation to:

export.tech@alifax.com

Alifax S.r.l.
Via Merano 30 - 33045 Nimis (UD) - Italy
Tel. +39-0432-197900
Fax: +39-0432-547378

IMPORTANT NOTES:

- Warranty expires 12 months after the installation date written on the present Installation Report duly filled in all its parts and returned to Alifax. If a copy of the report is not received back within one month, warranty will be invalidated.
- It is Your responsibility to verify that the device is properly installed and operates correctly and safely by filling this installation report and to keep such installation reports together with the related tests performed.

DISTRIBUTOR:	COUNTRY:
INSTALLATION SITE:	INSTALLATION DATE:
INSTRUMENT: ☐ TEST1 Model: _____	SERIAL NUMBER:
SOFTWARE VERSION (please see the print out):	

This document is divided in two parts: the installation report module (part one) containing installation qualification and operative qualification checklist and the installation procedure; the first one must be returned to manufacturer.

- Installation module as well as check list MUST be sent back to manufacturer by email

Please register at <http://www.alifax.com/> in order to obtain the latest version of the USER MANUAL and SERVICE MANUAL related to the software version installed on the instrument.

Фото 7 Инсталляционная форма



EU DECLARATION OF CONFORMITY

Code: ESR_DECO_SI1952XX_TEST1_3-0_EN

Date: 2022-05-25

The Manufacturer:

Manufacturer Name/Registered trade Name	ALIFAX S.R.L.
Address	Via Merano 30 33045 NIMIS (UD) ITALY

declares on his own responsibility that the following devices:

Product Name and Product Code	TEST1 THL	SI 195.210/THL
	TEST1 BCL	SI 195.220/BCL
	TEST1 SDL	SI 195.230/SDL
	TEST1 YDL	SI 195.240/YDL
	TEST1 MDL	SI 195.250/MDL
	TEST1 XDL	SI 195.260/XDL
Intended Purpose	TEST1 is an in vitro diagnostic automatic analyzer for the quantitative determination of Erythrocyte Sedimentation Rate (ESR) in human blood samples with EDTA from adult and pediatric patients with suspected inflammation. TEST1 provides results to inform clinical management for serious and non-serious conditions that require further diagnostic investigations and evaluation of the clinical status. The physician performs his/her expert evaluation on the information provided by the analyzer thanks to his/her professional knowledge, competence and skills as per the local law.	
Basic UDI-DI	805604014SIR20.195.804DA	
Risk Class	Class A, according to IVDR 746/2017.	

comply with the General Safety and Performance Requirements of the Regulation 746/2017 related to the In Vitro Diagnostic Medical Devices (IVDR).

We declare also that the devices have been designed and manufactured in conformity to the below:

Harmonised Standards	EN ISO 13485:2016 A11:2021: Medical devices - Quality management systems - Requirements for regulatory purposes EN ISO 14971:2019: Medical devices - Application of risk management to medical devices EN ISO 15223-1:2021: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied EN ISO 18113 -1:2011: In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) EN ISO 18113 -3:2011: Information supplied by the manufacturer (labelling) - Part 3: In vitro diagnostic instruments for professional use EN 62304:2006+A1:2015: Medical device software - Software life cycle processes
Technical Standards	IEC 62366-1:2015: Application of usability engineering to medical devices IEC TR 62366-2:2016: Guidance on the application of usability engineering to medical device IEC 61326-1:2012 Electrical equipment for measurement, control and laboratory use - EMC requirements IEC 61326-2-6:2012: Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 2-6: Particular requirements - In vitro diagnostic (IVD) medical equipment

DPGT0130 - Ed. 1 Rev. 1 - 2022-03-30

Pag. 1 di 2

Entered in Padova companies register at the n. 04337640280. Company stock 10.000.000 entirely deposited.

Company with single partner and subjected to direction and coordination of Alifax Holding S.p.A

©Copyright 2022. Documento segreto non divulgabile a terzi, di proprietà esclusiva di ALIFAX S.r.l. Ogni diritto riservato.

©Copyright 2022. This document is exclusive property of ALIFAX S.r.l, confidential and not freely disclosable. All rights reserved.

Фото 8 Декларация соответствия



FASCICOLO APPARECCHIO D195.05.03.02 Rev 1.9

Pag 8 di 9

DEVICE HISTORY RECORD

TEST1 THL / BCL / SDL / YDL / MDL / XDL / HDL (SI 195.2XX)

s.n. T17232commessa (Job Order) 2023-217200Data PRD 22/2/23

CERTIFICATO DI COLLAUDO TEST1		TEST1 - FINAL TEST REPORT	
MATRICOLA COMPONENTI - INTERNAL PARTS SERIAL NUMBERS			
DESCRIZIONE BOARD NAME	CODICE PART CODE	MATRICOLA SERIAL NUMBER	Sigla Operat. Operator
Lotto TELAIO	M195.11.13	-----124/2022-----	Ref
Lotto PANNELLO COMANDI	195.042	-----18/2022-----	
Lotto KIT FAST	195.032	-----01/2023-----	
Lotto ROTORE	195.031	-----103/2022-----	
Lotto TRASFORMATORE	195.046	-----1198/22-----	
SCHEDA ALIMENTATORE Power supply board	195.005	 195005*8865	Ref
SCHEDA CPU CPU Board	195.052	 195052*01324	Ref
SCHEDA CPS-MC CPS-MC Board	G1102.004	 G1102004*429	Ref
SCHEDA SENS. RACK E HOME Rack and Home Sensor Board	195.009	cos 38220071	Ref
SCHEDA SERIALE Serial Connectors Board	195.011	 195011*6066	Ref
CARRELLO Carriage	G195.001	 G195001*1675	Ref
LETTORE INTERNO DI CODICI A BARRE Internal barcode reader	ORTC1200- 1000	Serial Number G22M43586	Ref
COLLAUDI / FINAL TEST REPORT <ul style="list-style-type: none">TARATURA CPS-MC/ CPS-MC CalibrationPROVA DI FUNZIONAMENTO OPERATIVO / Final TestINSERIMENTO E VALUTAZIONE DEI CORRETTI PARAMETRI DI SETUP / Setup default parameters check <p>SI CERTIFICA CHE TUTTI I CONTROLLI INTERNI SONO STATI SUPERATI CON ESITO POSITIVO WE CERTIFY THAT ALL INTERNAL CHECKS WERE POSITIVELY PASSED</p> <p>Approvato dal Responsabile Produzione Strumenti / Manufacturer Manager Approval</p> <p>DATA/Date <u>16/5/2023</u> FIRMA/Signature </p>			

©Copyright 2022. Documento segreto non divulgabile a terzi, di proprietà esclusiva di ALIFAX S.r.l. Ogni diritto riservato.
©Copyright 2022. This document is exclusive property of ALIFAX S.r.l, confidential and not freely disclosable. All rights reserved.

Фото 9 Сертификат выходного контроля качества (анализатор TEST1 модель THL)



Фото 10 Подставка под штативы





Фото 12 Смарт-карты TEST1 на 1000/4000/10000/20000 тестов




Фото 13 Штатив металлический Alifax (только для TEST1 THL)




Фото 14 Внешний сканер штрих-кодов (только для TEST1 THL)



Фото 15 Крышка-адаптер для пробирок Sarstedt Monovette (только для TEST1 THL)



Фото 16 Предохранители



Фото 17 Кабель питания



Фото 18 Стартовый набор термобумаги в рулонах



Wash tank
Réservoir de lavage

Min level 100ml
It is recommended to fill
the wash tank before reaching the
black line level

Niveau minimum 100 ml
On recommande de remplir avec
de l'eau distillée avant de rejoindre
le niveau indiqué par la ligne

Фото 19 Резервуар для промывки



Фото 20 Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 BCL (общий вид)



Фото 21 Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 BCL (фото маркировки)



TEST 1

Руководство по эксплуатации

Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 с принадлежностями

Действительно для всех моделей

Применимо к версии 9.00X



Редакция для Российской Федерации



**Прибор для профессионального использования в области
In Vitro диагностики**

Copyright © Alifax S.r.l.

Информация в данном руководстве является собственностью ALIFAX. Все права защищены.

Несанкционированное копирование руководства или его частей запрещено.

Торговая марка TEST 1 является собственностью Alifax S.r.l.



Фото 23 Руководство по эксплуатации на электронном носителе (CD-диск)



Фото 24 Руководство по эксплуатации на электронном носителе (флеш-карта)

FASCICOLO APPARECCHIO D195.05.03.02 Rev 1.9
DEVICE HISTORY RECORD

TEST1 THL / BCL / SDL / YDL / MDL / XDL / HDL (SI 195.2XX)
 s.n. 777236 commessa (Job Order) 2023-2/1532

Pag. 9 di 9

Data PRD 28/12/23

CORREDO DI SERIE / PACKING LIST

STRUMENTO INSTRUMENT	MODELLO MODEL						LINGUA LANGUAGE	
	THL	BCL	SDL	YDL	MDL	XDL	EN	IT
☒ TEST1 SI 195.2XX	☐	☒	☐	☐	☐	☐	☒	☐

Corredo di serie / Packing list

- ☒ N.° 2 SI195800 Rotoli Carta termica / Thermal Paper rolls
- ☒ N.° 2 SI205801 Tanica di lavaggio / Wash Tank
- ☒ N.° 2 SI205801 Tanica di scarico / Waste Tank
- ☒ N.° 2 Fusibili di alimentazione T5A / Main fuses T5A
- ☒ N.° 1 Cavo Alimentazione / Main cable plug
- ☐ Cavo alimentazione personalizzato/Personalized main cable _____
- ☒ N.° 1 Ago a vite / Needle
- ☒ N.° 1 CD Manuale operativo multilingua / CD Multi-Language Operative manual (Language and version)
 ESR_UM_SI195_TEST1_SW 11-01 A-A-A-EN-ES-IT-PT-FR-DE-FI-LT ESR-CD-SI195-TEST1-SW 11-00X-1-2
- ☒ ESR_IQOQ_SI195_TEST1: Modulo di installazione / Installation Report

SOLO STRUMENTI SDL / SDL INSTRUMENTS ONLY:

- ☐ N.° 4 Pettine rack Sysmex XN / Sysmex XN rack scallop
- ☐ Nota tecnica n.° 2004 EN / Technical note n.° 2004 EN

SOLO STRUMENTI SDL PER PROVETTE SARSTEDT / SDL INSTRUMENTS FOR SARSTEDT TUBES ONLY:

- ☐ N.° 4 Pettine x Sarstedt 11X66 / Sarstedt 11x66 scallop
- ☐ Nota tecnica n.° 2006 EN / Technical note n.° 2006 EN

- ☒ Dichiarazione di conformità CE / EC declaration of conformity
- ☒ DPRD0105: Prove elettriche di routine / Electrical routine tests
- ☒ Pag 8/9: Certificato di Collaudo / Final Test report
- ☒ Pag 9/9: Packing List / This document

Compilatore/Signature Cossiga

Data/Date 21/3/2023

Approvato dal Responsabile Produzione Strumenti per rilascio finale / Manufacturer Manager Approval

DATA/Date 24/03/2023 FIRMA/Signature _____

Approvato da Assicurazione Qualità (Approved by Quality Assurance)

Data 21/03/2023 Firma _____
 (Date) (Signature)

ALIFAX S.r.l.
 via Petrarca 2/1 – 35020 POLVERARA (PD) – ITALY
 Tel. +39 0490992000 e-mail: info@alifax.com
 PI 04337640280



ALIFAX S.r.l.
 via Merano 30 – 33045 NIMIS (UD) – ITALY
 Documents available on: www.alifax.com

©Copyright 2022. Documento segreto non divulgabile a terzi, di proprietà esclusiva di ALIFAX S.r.l. Ogni diritto riservato.
 ©Copyright 2022. This document is exclusive property of ALIFAX S.r.l, confidential and not freely disclosable. All rights reserved.

Фото 25 Упаковочный лист (анализатор TEST1 модель BCL)



INSTALLATION REPORT for TEST1 Instruments

Rev. 3.0 - Valid since 2022 May 25

Dear Distributor,

You are kindly requested to fill-in this form during the installation and return a copy within one month since the installation to:

export.tech@alifax.com

Alifax S.r.l.
Via Merano 30 - 33045 Nimis (UD) - Italy
Tel. +39-0432-197900
Fax: +39-0432-547378

IMPORTANT NOTES:

- Warranty expires 12 months after the installation date written on the present Installation Report duly filled in all its parts and returned to Alifax. If a copy of the report is not received back within one month, warranty will be invalidated.
- It is Your responsibility to verify that the device is properly installed and operates correctly and safely by filling this installation report and to keep such installation reports together with the related tests performed.

DISTRIBUTOR:	COUNTRY:
INSTALLATION SITE:	INSTALLATION DATE:
INSTRUMENT: ☐ TEST1 Model: _____	SERIAL NUMBER:
SOFTWARE VERSION (please see the print out):	

This document is divided in two parts: the installation report module (part one) containing **installation qualification** and **operative qualification** checklist and the installation procedure; the first one must be returned to manufacturer.

- Installation module as well as check list MUST be sent back to manufacturer by email

Please register at <http://www.alifax.com/> in order to obtain the latest version of the USER MANUAL and SERVICE MANUAL related to the software version installed on the instrument.

Фото 26 Инсталляционная форма



EU DECLARATION OF CONFORMITY

Code: ESR_DECO_SI1952XX_TEST1_3-0_EN

Date: 2022-05-25

The Manufacturer:

Manufacturer Name/Registered trade Name	ALIFAX S.R.L.
Address	Via Merano 30 33045 NIMIS (UD) ITALY

declares on his own responsibility that the following devices:

Product Name and Product Code	TEST1 THL	SI 195.210/THL
	TEST1 BCL	SI 195.220/BCL
	TEST1 SDL	SI 195.230/SDL
	TEST1 YDL	SI 195.240/YDL
	TEST1 MDL	SI 195.250/MDL
	TEST1 XDL	SI 195.260/XDL
Intended Purpose	TEST1 is an in vitro diagnostic automatic analyzer for the quantitative determination of Erythrocyte Sedimentation Rate (ESR) in human blood samples with EDTA from adult and pediatric patients with suspected inflammation. TEST1 provides results to inform clinical management for serious and non-serious conditions that require further diagnostic investigations and evaluation of the clinical status. The physician performs his/her expert evaluation on the information provided by the analyzer thanks to his/her professional knowledge, competence and skills as per the local law.	
Basic UDI-DI	805604014SIR20.195.804DA	
Risk Class	Class A, according to IVDR 746/2017.	

comply with the General Safety and Performance Requirements of the Regulation 746/2017 related to the In Vitro Diagnostic Medical Devices (IVDR).

We declare also that the devices have been designed and manufactured in conformity to the below:

Harmonised Standards	EN ISO 13485:2016 A11:2021: Medical devices - Quality management systems - Requirements for regulatory purposes EN ISO 14971:2019: Medical devices - Application of risk management to medical devices EN ISO 15223-1:2021: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied EN ISO 18113 -1:2011: In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) EN ISO 18113 -3:2011: Information supplied by the manufacturer (labelling) - Part 3: In vitro diagnostic instruments for professional use EN 62304:2006+A1:2015: Medical device software - Software life cycle processes
Technical Standards	IEC 62366-1:2015: Application of usability engineering to medical devices IEC TR 62366-2:2016: Guidance on the application of usability engineering to medical device IEC 61326-1:2012 Electrical equipment for measurement, control and laboratory use - EMC requirements IEC 61326-2-6:2012: Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 2-6: Particular requirements - In vitro diagnostic (IVD) medical equipment

DPGT0130 – Ed. 1 Rev. 1 – 2022-03-30

Pag. 1 di 2

Entered in Padova companies register at the n. 04337640280. Company stock 10.000.000 entirely deposited.
Company with single partner and subjected to direction and coordination of Alifax Holding S.p.A.
©Copyright 2022. Documento segreto non divulgabile a terzi, di proprietà esclusiva di ALIFAX S.r.l. Ogni diritto riservato.
©Copyright 2022. This document is exclusive property of ALIFAX S.r.l, confidential and not freely disclosable. All rights reserved.

Фото 27 Декларация соответствия

CERTIFICATO DI COLLAUDO TEST1

TEST1 - FINAL TEST REPORT

MATRICOLA COMPONENTI - INTERNAL PARTS SERIAL NUMBERS

DESCRIZIONE BOARD NAME	CODICE PART CODE	MATRICOLA SERIAL NUMBER	Sigla Operat. Operator
Lotto TELAIO	M195.11.13	<u>124/2022</u>	
Lotto PANNELLO COMANDI	195.042	<u>02/2023</u>	
Lotto KIT FAST	195.032	<u>02/2023</u>	
Lotto ROTORE	195.031	<u>103/2022</u>	
Lotto TRASFORMATORE	195.046	<u>1138/22</u>	
SCHEDA ALIMENTATORE Power supply board	195.005	 195005*8942	
SCHEDA CPU CPU Board	195.052	 195052*01423	
SCHEDA CPS-MC CPS-MC Board	G1102.004	 G1102004*457	
SCHEDA SENS. RACK E HOME Rack and Home Sensor Board	195.009	<u>00903230177</u>	
SCHEDA SERIALE Serial Connectors Board	195.011	 195011*6111	
CARRELLO Carriage	G195.001	 G195001*1687	
LETTORE INTERNO DI CODICI A BARRE Internal barcode reader	ORTC1200- 1000	 Serial Number G22M43556	

COLLAUDI / FINAL TEST REPORT

- TARATURA CPS-MC/ CPS-MC Calibration
- PROVA DI FUNZIONAMENTO OPERATIVO / Final Test
- INSERIMENTO E VALUTAZIONE DEI CORRETTI PARAMETRI DI SETUP / Setup default parameters check

SI CERTIFICA CHE TUTTI I CONTROLLI INTERNI SONO STATI SUPERATI CON ESITO POSITIVO
WE CERTIFY THAT ALL INTERNAL CHECKS WERE POSITIVELY PASSED

Approvato dal Responsabile Produzione Strumenti / Manufacturer Manager Approval


DATA/Date 21/05/2023 FIRMA/Signature 



Фото 29 Подставка под штативы





Фото 31 Смарт-карты TEST1 на 1000/4000/10000/20000 тестов




Фото 32 Штатив Alifax для TEST1 BCL (зеленый)



Фото 33 Кабель питания



Фото 34 Стартовый набор термобумаги в рулонах



Фото 35 Предохранители



Фото 36 Резервуар для промывки



Фото 37 Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 SDL (общий вид)



Фото 38 Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 SDL (фото маркировки)



TEST 1

Руководство по эксплуатации


Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 с принадлежностями

Действительно для всех моделей

Применимо к версии 9.00X



Редакция для Российской Федерации



**Прибор для профессионального использования в области
In Vitro диагностики**

Copyright © Alifax S.r.l.

Информация в данном руководстве является собственностью ALIFAX. Все права защищены.
Несанкционированное копирование руководства или его частей запрещено.
Торговая марка TEST 1 является собственностью Alifax S.r.l.



Фото 40 Руководство по эксплуатации на электронном носителе (CD-диск)



Фото 41 Руководство по эксплуатации на электронном носителе (флеш-карта)



FASCICOLO APPARECCHIO D195.05.03.02 Rev 1.9

Pag. 2 di 2

DEVICE HISTORY RECORD

TEST1 THL / BCL / SDL / YDL / MDL / XDL / HDL (SI 195.2XX)

s.n. 717223commessa (Job Order) 2023 2/000Data PRD 16/12/2023

CORREDO DI SERIE / PACKING LIST

STRUMENTO INSTRUMENT	MODELLO MODEL						LINGUA LANGUAGE	
	THL	BCL	SDL	YDL	MDL	XDL	EN	IT
☒ TEST1 SI 195.2XX	☐	☐	☒	☐	☐	☐	☒	☐

Corredo di serie / Packing list

- ☒ N.° 2 SI195800 Rotoli Carta termica / Thermal Paper rolls
- ☒ N.° 2 SI205801 Tanica di lavaggio / Wash Tank
- ☒ N.° 2 SI205801 Tanica di scarico / Waste Tank
- ☒ N.° 2 Fusibili di alimentazione T5A / Main fuses T5A
- ☒ N.° 1 Cavo Alimentazione / Main cable plug
- ☐ Cavo alimentazione personalizzato/Personalized main cable _____
- ☒ N.° 1 Ago a vite / Needle
- ☒ N.° 1 CD Manuale operativo multilingua / CD Multi-Language Operative manual (Language and version)
ESR_UM_SI195_TEST1_SW 11-00A-1-0-EN-ES-IT-PT-FR-DE-FI ESR-CD-SI195-TEST1-SW 11-COX-1-1
- ☒ ESR_IQOQ_SI195_TEST1: Modulo di installazione / Installation Report

SOLO STRUMENTI SDL / SDL INSTRUMENTS ONLY:

- ☒ N.° 4 Pettine rack Sysmex XN / Sysmex XN rack scallop
- ☒ Nota tecnica n.° 2004 EN / Technical note n.° 2004 EN

SOLO STRUMENTI SDL PER PROVETTE SARSTEDT / SDL INSTRUMENTS FOR SARSTEDT TUBES ONLY:

- ☐ N.° 4 Pettine x Sarstedt 11X66 / Sarstedt 11x66 scallop
- ☐ Nota tecnica n.° 2006 EN / Technical note n.° 2006 EN

- ☒ Dichiarazione di conformità CE / EC declaration of conformity
- ☒ DPRD0105: Prove elettriche di routine / Electrical routine tests
- ☒ Pag 8/9: Certificato di Collaudo / Final Test report
- ☒ Pag 9/9: Packing List / This document

Compilatore/Signature CossigaData/Date 08/03/2023

Approvato dal Responsabile Produzione Strumenti per rilascio finale / Manufacturer Manager Approval

DATA/Date 08/03/2023 FIRMA/Signature _____

Approvato da Assicurazione Qualità (Approved by Quality Assurance)

Data 08/03/2023 Firma _____
(Date)

(Signature)

ALIFAX S.r.l.
via Petrarca 2/1 - 35020 POLVERARA (PD) - ITALY
Tel. +39 0490992000 e-mail: info@alifax.com
PI 04337640280



ALIFAX S.r.l.
via Merano 30 - 33045 NIMIS (UD) - ITALY
Documents available on: www.alifax.com

©Copyright 2022. Documento segreto non divulgabile a terzi, di proprietà esclusiva di ALIFAX S.r.l. Ogni diritto riservato.
©Copyright 2022. This document is exclusive property of ALIFAX S.r.l, confidential and not freely disclosable. All rights reserved.

Фото 42 Упаковочный лист (анализатор TEST1 модель SDL)



INSTALLATION REPORT for TEST1 Instruments

Rev. 3.0 - Valid since 2022 May 25

Dear Distributor,

You are kindly requested to fill-in this form during the installation and return a copy within one month since the installation to:

export.tech@alifax.com

Alifax S.r.l.
Via Merano 30 - 33045 Nimis (UD) - Italy
Tel. +39-0432-197900
Fax: +39-0432-547378

IMPORTANT NOTES:

- Warranty expires 12 months after the installation date written on the present Installation Report duly filled in all its parts and returned to Alifax. If a copy of the report is not received back within one month, warranty will be invalidated.
- It is Your responsibility to verify that the device is properly installed and operates correctly and safely by filling this installation report and to keep such installation reports together with the related tests performed.


DISTRIBUTOR:	COUNTRY:
INSTALLATION SITE:	INSTALLATION DATE:
INSTRUMENT: ☐ TEST1 Model: _____	SERIAL NUMBER:
SOFTWARE VERSION (please see the print out):	

This document is divided in two parts: the installation report module (part one) containing **installation qualification** and **operative qualification** checklist and the installation procedure; the first one must be returned to manufacturer.

- Installation module as well as check list MUST be sent back to manufacturer by email

Please register at <http://www.alifax.com/> in order to obtain the latest version of the USER MANUAL and SERVICE MANUAL related to the software version installed on the instrument.

Фото 43 Инсталляционная форма



EU DECLARATION OF CONFORMITY

Code: ESR_DECO_SI1952XX_TEST1_3-0_EN

Date: 2022-05-25

The Manufacturer:

Manufacturer Name/Registered trade Name	ALIFAX S.R.L.
Address	Via Merano 30 33045 NIMIS (UD) ITALY

declares on his own responsibility that the following devices:

Product Name and Product Code	TEST1 THL	SI 195.210/THL
	TEST1 BCL	SI 195.220/BCL
	TEST1 SDL	SI 195.230/SDL
	TEST1 YDL	SI 195.240/YDL
	TEST1 MDL	SI 195.250/MDL
	TEST1 XDL	SI 195.260/XDL
Intended Purpose	TEST1 is an in vitro diagnostic automatic analyzer for the quantitative determination of Erythrocyte Sedimentation Rate (ESR) in human blood samples with EDTA from adult and pediatric patients with suspected inflammation. TEST1 provides results to inform clinical management for serious and non-serious conditions that require further diagnostic investigations and evaluation of the clinical status. The physician performs his/her expert evaluation on the information provided by the analyzer thanks to his/her professional knowledge, competence and skills as per the local law.	
Basic UDI-DI	805604014SIR20.195.804DA	
Risk Class	Class A, according to IVDR 746/2017.	

comply with the General Safety and Performance Requirements of the Regulation 746/2017 related to the In Vitro Diagnostic Medical Devices (IVDR).

We declare also that the devices have been designed and manufactured in conformity to the below:

Harmonised Standards	EN ISO 13485:2016 A11:2021: Medical devices - Quality management systems - Requirements for regulatory purposes EN ISO 14971:2019: Medical devices - Application of risk management to medical devices EN ISO 15223-1:2021: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied EN ISO 18113 -1:2011: In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) EN ISO 18113 -3:2011: Information supplied by the manufacturer (labelling) - Part 3: In vitro diagnostic instruments for professional use EN 62304:2006+A1:2015: Medical device software - Software life cycle processes
Technical Standards	IEC 62366-1:2015: Application of usability engineering to medical devices IEC TR 62366-2:2016: Guidance on the application of usability engineering to medical device IEC 61326-1:2012 Electrical equipment for measurement, control and laboratory use - EMC requirements IEC 61326-2-6:2012: Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 2-6: Particular requirements - In vitro diagnostic (IVD) medical equipment

DPGT0130 - Ed. 1 Rev. 1 - 2022-03-30

Pag. 1 di 2

Entered in Padova companies register at the n. 04337640280. Company stock 10.000.000 entirely deposited.

Company with single partner and subjected to direction and coordination of Alifax Holding S.p.A.

©Copyright 2022. Documento segreto non divulgabile a terzi, di proprietà esclusiva di ALIFAX S.r.l. Ogni diritto riservato.

©Copyright 2022. This document is exclusive property of ALIFAX S.r.l, confidential and not freely disclosable. All rights reserved.

Фото 44 Декларация соответствия



FASCICOLO APPARECCHIO D196.05.03.02 Rev 1.9

Pag. 8/8

DEVICE HISTORY RECORD

TEST1 THL / BCL / SDL / YDL / MDL / XDL / HDL (SI 195.2XX)

s.n. 717220 commessa (Job Order) 2023-21930Data PRD 16/12/2023

CERTIFICATO DI COLLAUDO TEST1		TEST1 - FINAL TEST REPORT	
MATRICOLA COMPONENTI - INTERNAL PARTS SERIAL NUMBERS			
DESCRIZIONE BOARD NAME	CODICE PART CODE	MATRICOLA SERIAL NUMBER	Sigla Operat. Operator
Lotto TELAIO	M195.11.13	<u>12/2022</u>	
Lotto PANNELLO COMANDI	195.042	<u>12/2022</u>	
Lotto KIT FAST	195.032	<u>20/2022</u>	
Lotto ROTORE	195.031	<u>10/2022</u>	
Lotto TRASFORMATORE	195.046	<u>19/8/22</u>	
SCHEDA ALIMENTATORE Power supply board	195.005	 195005*8869	
SCHEDA CPU CPU Board	195.052	 195052*01218	
SCHEDA CPS-MC CPS-MC Board	G1102.004	 G1102004*406	
SCHEDA SENS. RACK E HOME Rack and Home Sensor Board	195.009	<u>00903230161</u>	
SCHEDA SERIALE Serial Connectors Board	195.011	 195011*5744	
CARRELLO Carriage	G195.001	 G195001*1680	
LETTORE INTERNO DI CODICI A BARRE Internal barcode reader	ORTC1200- 1000	Serial Number G22M43598	
COLLAUDI / FINAL TEST REPORT <ul style="list-style-type: none">TARATURA CPS-MC/ CPS-MC CalibrationPROVA DI FUNZIONAMENTO OPERATIVO / Final TestINSERIMENTO E VALUTAZIONE DEI CORRETTI PARAMETRI DI SETUP / Setup default parameters check <p>SI CERTIFICA CHE TUTTI I CONTROLLI INTERNI SONO STATI SUPERATI CON ESITO POSITIVO WE CERTIFY THAT ALL INTERNAL CHECKS WERE POSITIVELY PASSED</p> <p>Approvato dal Responsabile Produzione Strumenti / Manufacturer Manager Approval</p> <p>DATA/Date <u>08/01/2023</u> FIRMA/Signature </p>			



Фото 46 Подставка под штативы





Фото 48 Смарт-карты TEST1 на 1000/4000/10000/20000 тестов




Фото 49 Штатив Alifax для TEST1 SDL (желтый)



Фото 50 Кабель питания



Фото 51 Стартовый набор термобумаги в рулонах



Фото 52 Предохранители



Фото 53 Резервуар для промывки



Фото 54 Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 YDL (общий вид)



Фото 55 Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 YDL (фото маркировки)



TEST 1

Руководство по эксплуатации

Анализатор автоматический для определения скорости оседания эритроцитов (СОЭ) TEST1 с принадлежностями

Действительно для всех моделей

Применимо к версии 9.00X



Редакция для Российской Федерации



**Прибор для профессионального использования в области
In Vitro диагностики**

Copyright © Alifax S.r.l.

Информация в данном руководстве является собственностью ALIFAX. Все права защищены.
Несанкционированное копирование руководства или его частей запрещено.
Торговая марка TEST 1 является собственностью Alifax S.r.l.



Фото 57 Руководство по эксплуатации на электронном носителе (CD-диск)



Фото 58 Руководство по эксплуатации на электронном носителе (флеш-карта)



FASCICOLO APPARECCHIO D195.05.03.02 Rev 1.7

DEVICE HISTORY RECORD

Pag. 9 di 9

TEST1 THL / BCL / SDL / YDL / MDL / XDL / HDL (SI 195.2XX)

s.n. T46865

commessa (Job Order) 2022-2/773

Data PRD 14/02/22

CORREDO DI SERIE / PACKING LIST

STRUMENTO INSTRUMENT	MODELLO MODEL						LINGUA LANGUAGE	
	THL	BCL	SDL	YDL	MDL	XDL	EN	IT
☒ TEST1 SI 195.2XX	☐	☐	☐	☒	☐	☐	☒	☐

Corredo di serie / Packing list

- ☒ N.° 2 SI195800 Rotoli Carta termica / Thermal Paper rolls
- ☒ N.° 2 SI205801 Tanica di lavaggio / Wash Tank
- ☒ N.° 2 SI205801 Tanica di scarico / Waste Tank
- ☒ N.° 2 Fusibili di alimentazione T5A / Main fuses T5A
- ☒ N.° 1 Cavo Alimentazione / Main cable plug
- ☐ Cavo alimentazione personalizzato/Personalized main cable _____
- ☒ N.° 1 Ago a vite / Needle
- ☒ N.° 1 CD Manuale operativo multilingua / CD Multi-Language Operative manual (Language and version)
ESR_UM_SI195_TEST1_SW 9-00C-A-9-EN-ES-IT-PT
- ☒ ESR_IQOQ_SI195_TEST1: Modulo di installazione / Installation Report

SOLO STRUMENTI SDL / SDL INSTRUMENTS ONLY:

- ☐ N.° 4 Pettine rack Sysmex XN / Sysmex XN rack scallop
- ☐ Nota tecnica n.° 2004 EN / Technical note n.° 2004 EN

SOLO STRUMENTI SDL PER PROVETTE SARSTEDT / SDL INSTRUMENTS FOR SARSTEDT TUBES ONLY:

- ☐ N.° 4 Pettine x Sarstedt 11X66 / Sarstedt 11x66 scallop
- ☐ Nota tecnica n.° 2006 EN / Technical note n.° 2006 EN

- ☒ Dichiarazione di conformità CE / EC declaration of conformity
- ☒ DPRD0105: Prove Elettriche EN61010-2-101:2002 Allegato F / Electrical routine tests EN61010-2-101:2002 Annex F
- ☒ Pag 8/9: Certificato di Collaudo / Final Test report
- ☒ Pag 9/9: Packing List / This document

Compilatore/Signature CossichData/Date 30/03/2022

Approvato dal Responsabile Produzione Strumenti per rilascio finale / Manufacturer Manager Approval

DATA/Date 30/03/2022FIRMA/Signature [Signature]

Approvato da Assicurazione Qualità (Approved by Quality Assurance)

Data 30/03/2022
(Date)Firma [Signature]
(Signature)

ALIFAX S.r.l.
via Petrarca 2/1 - 35020 POLVERARA (PD) - ITALY
Tel. +39 0490992000 e-mail: info@alifax.com
PI 04337640280



ALIFAX S.r.l.
via Merano 30 - 33045 NIMIS (UD) - ITALY
Documents available on: www.alifax.com

©Copyright 2022. Documento segreto non divulgabile a terzi, di proprietà esclusiva di ALIFAX S.r.l. Ogni diritto riservato.
©Copyright 2022. This document is exclusive property of ALIFAX S.r.l, confidential and not freely disclosable. All rights reserved.

Фото 59 Упаковочный лист (анализатор TEST1 модель YDL)



INSTALLATION REPORT for TEST1 Instruments

Rev. 3.0 - Valid since 2022 May 25

Dear Distributor,

You are kindly requested to fill-in this form during the installation and return a copy within one month since the installation to:

export.tech@alifax.com

Alifax S.r.l.
Via Merano 30 - 33045 Nimis (UD) - Italy
Tel. +39-0432-197900
Fax: +39-0432-547378

IMPORTANT NOTES:

- Warranty expires 12 months after the installation date written on the present Installation Report duly filled in all its parts and returned to Alifax. If a copy of the report is not received back within one month, warranty will be invalidated.
- It is Your responsibility to verify that the device is properly installed and operates correctly and safely by filling this installation report and to keep such installation reports together with the related tests performed.


DISTRIBUTOR:	COUNTRY:
INSTALLATION SITE:	INSTALLATION DATE:
INSTRUMENT: ☐ TEST1 Model: _____	SERIAL NUMBER:
SOFTWARE VERSION (please see the print out):	

This document is divided in two parts: the installation report module (part one) containing **installation qualification** and **operative qualification** checklist and the installation procedure; the first one must be returned to manufacturer.

- Installation module as well as check list MUST be sent back to manufacturer by email

Please register at <http://www.alifax.com/> in order to obtain the latest version of the USER MANUAL and SERVICE MANUAL related to the software version installed on the instrument.

Фото 60 Инсталляционная форма



EU DECLARATION OF CONFORMITY

Code: ESR_DECO_SI1952XX_TEST1_3-0_EN

Date: 2022-05-25

The Manufacturer:

Manufacturer Name/Registered trade Name	ALIFAX S.R.L.
Address	Via Merano 30 33045 NIMIS (UD) ITALY

declares on his own responsibility that the following devices:

Product Name and Product Code	TEST1 THL	SI 195.210/THL
	TEST1 BCL	SI 195.220/BCL
	TEST1 SDL	SI 195.230/SDL
	TEST1 YDL	SI 195.240/YDL
	TEST1 MDL	SI 195.250/MDL
	TEST1 XDL	SI 195.260/XDL
Intended Purpose	TEST1 is an in vitro diagnostic automatic analyzer for the quantitative determination of Erythrocyte Sedimentation Rate (ESR) in human blood samples with EDTA from adult and pediatric patients with suspected inflammation. TEST1 provides results to inform clinical management for serious and non-serious conditions that require further diagnostic investigations and evaluation of the clinical status. The physician performs his/her expert evaluation on the information provided by the analyzer thanks to his/her professional knowledge, competence and skills as per the local law.	
Basic UDI-DI	805604014SIR20.195.804DA	
Risk Class	Class A, according to IVDR 746/2017.	

comply with the General Safety and Performance Requirements of the Regulation 746/2017 related to the In Vitro Diagnostic Medical Devices (IVDR).

We declare also that the devices have been designed and manufactured in conformity to the below:

Harmonised Standards	EN ISO 13485:2016 A11:2021: Medical devices - Quality management systems - Requirements for regulatory purposes EN ISO 14971:2019: Medical devices - Application of risk management to medical devices EN ISO 15223-1:2021: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied EN ISO 18113 -1:2011: In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) EN ISO 18113 -3:2011: Information supplied by the manufacturer (labelling) - Part 3: In vitro diagnostic instruments for professional use EN 62304:2006+A1:2015: Medical device software - Software life cycle processes
Technical Standards	IEC 62366-1:2015: Application of usability engineering to medical devices IEC TR 62366-2:2016: Guidance on the application of usability engineering to medical device IEC 61326-1:2012 Electrical equipment for measurement, control and laboratory use - EMC requirements IEC 61326-2-6:2012: Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 2-6: Particular requirements - In vitro diagnostic (IVD) medical equipment

DPGT0130 - Ed. 1 Rev. 1 - 2022-03-30

Pag. 1 di 2

Entered in Padova companies register at the n. 04337640280. Company stock 10.000.000 entirely deposited.
Company with single partner and subjected to direction and coordination of Alifax Holding S.p.A.
©Copyright 2022. Documento segreto non divulgabile a terzi, di proprietà esclusiva di ALIFAX S.r.l. Ogni diritto riservato.
©Copyright 2022. This document is exclusive property of ALIFAX S.r.l, confidential and not freely disclosable. All rights reserved.

Фото 61 Декларация соответствия



CERTIFICATO DI COLLAUDO TEST1		TEST1 - FINAL TEST REPORT	
MATRICOLA COMPONENTI - INTERNAL PARTS SERIAL NUMBERS			
DESCRIZIONE BOARD NAME	CODICE PART CODE	MATRICOLA SERIAL NUMBER	Sigla Operat. Operator
Lotto TELAIO	M195.11.13	83/22	2
Lotto PANNELLO COMANDI	195.042	01/22	
Lotto KIT FAST	195.032	01/22	
Lotto ROTORE	195.031	14/22	
Lotto TRASFORMATORE	195.046	1299/20	
SCHEDA ALIMENTATORE Power supply board	195.005	 195005*8273	2
SCHEDA CPU CPU Board	195.052	 195052*01401	2
SCHEDA CPS-MC CPS-MC Board	G1102.003	 G1102003*1287	2
SCHEDA SENS. RACK E HOME Rack and Home Sensor Board	195.009	00926210118	2
SCHEDA SERIALE Serial Connectors Board	195.011	 195011*5558	2
CARRELLO Carriage	G195.001	 G195001*1310	2
LETTORE INTERNO DI CODICI A BARRE Internal barcode reader	ORTC1200- 1000	Serial Number G21M40675	2
COLLAUDI / FINAL TEST REPORT <ul style="list-style-type: none">TARATURA CPS-MC/ CPS-MC CalibrationPROVA DI FUNZIONAMENTO OPERATIVO / Final TestINSERIMENTO E VALUTAZIONE DEI CORRETTI PARAMETRI DI SETUP / Setup default parameters check <p>SI CERTIFICA CHE TUTTI I CONTROLLI INTERNI SONO STATI SUPERATI CON ESITO POSITIVO WE CERTIFY THAT ALL INTERNAL CHECKS WERE POSITIVELY PASSED</p> <p>Approvato dal Responsabile Produzione Strumenti / Manufacturer Manager Approval</p> <p>DATA/Date <u>30/03/2022</u> FIRMA/Signature <u>[Signature]</u></p>			



Фото 63 Подставка под штативы





Фото 65 Смарт-карты TEST1 на 1000/4000/10000/20000 тестов




Фото 66 Штатив Alifax для TEST1 YDL (синий)



Фото 67 Кабель питания



Фото 68 Стартовый набор термобумаги в рулонах



Фото 69 Предохранители



Wash tank
Réservoir de lavage

Min level 100ml
It is recommended to fill
the wash tank before reaching the
black line level

Niveau minimum 100 ml
On recommande de remplir avec
de l'eau distillée avant de rejoindre
le niveau indiqué par la ligne

В данном документе
пронумеровано и прошито

71

лист(а)ов

