


REF Q-NCOV-01C

STANDARD Q

COVID-19 IgM/IgG Combo

CONTENTS

- 1) Test device (individually in a foil pouch with desiccant) x 40
- 2) Buffer bottle x 1
- 3) Capillary tube (20µl) x 40
- 4) Instructions for use x 1



www.sdbiosensor.com

Q.C

 **SD BIOSENSOR**

STANDARD Q
COVID-19 IgM/IgG

CONTENTS

- 1) Test device (prepack)
- 2) Buffer solution
- 3) Capillary tube (pipette)
- 4) Instructions for use



www.goszdraznadzor.ru

STANDARD Q

COVID-19 IgM/IgG Combo Test



LOT NO.:



MFG DATE:



EXP DATE:

QCO5020011P / Sub : P-1

-Buffer : COD220026


2020.07.06.

2022.01.05.

RECOMMENDED TECHNIQUES TO IMPROVE PERFORMANCE

Technique 1
Drop the buffer

Technique 2.
Attach the provided film on to the test device.



STANDARD Q
Rapid Test

CE IVD

STANDARD

Capillary tube (20µl)

LOT No.: BS202006
EXP: 2025.06.10,
Quantity: 20PCS

SD BIOSENSOR


LOT No.: BS202006
EXP: 2025.06.10
Quantity: 20PCS

SD BIOSENSOR



PI Please attach
P performing
C circular

STANDARD Q
Rapid Test



STANDARDTM

/ Item **STANDARD Q COVID-19 IgM/IgG Combo**

REF / REF No. **Q-NCOV-01C**

LOT / LOT No. **QCO5020011P**

MFG DATE **2020.07.06.**

EXP DATE **2022.01.05.**

#Buffer LOT : COD220026

PL23COVSENRO
Issue date: 2020.04

Manufactured by



SD BIOSENSOR

에스디 바이오센서(주)
www.sdbiosensor.com

STANDARD Q COVID
Combo

Buffer bottle



COD220026



2020.05.29.



SD BIOSENSOR

Capillary tube (20μl)

LOT No. : BS202006

EXP : 2025.06.10.

Quantity : 20PCS



IVD



 **SD BIOSENSOR**

L46RT4ENR0

Issue date: 2020.04

Please attach this film on the result window of device, when performing the test outdoor or in an environment with lot of air circulation (clean bench).

Please attach this film on the result window of device, when performing the test outdoor or in an environment with lot of air circulation (clean bench).

L24RT2ENR0 / Issue date: 2020.04

Прошито, пронумеровано и
скреплено печатью 7 листов.

Генеральный директор
ООО «Уайт Продакт»
Немцов В.Б.



Макет этикетки потребительской упаковки

	Набор реагентов для иммунохроматографического выявления антител IgM/IgG к вирусу SARS-CoV-2 в образцах цельной крови, сыворотки или плазмы (STANDARD Q COVID-19 IgM/IgG Combo). Lot №QCO5020011P SD Biosensor, Inc. (ЭС-ДИ Биосенсор, Инк.) <i>Главный офис:</i> Республика Корея, 16690, провинция Кёнги-до, г. Сувон, Ёнтхон-гу, Тогён-дэро, 1556-киль, 16, этаж С-4~5 <i>Адрес производства:</i> Республика Корея, 28161, провинция Чхунчхон-пукто, г. Чхончжу, Хындок-гу, Осон-ып, Осон сэнмён-но-4, 74 Представитель производителя в РФ: ООО «Уайт Продакт»: 117105, г. Москва, 1-й Нагатинский пр-д, д. 2, стр. 7, этаж 3, помещение VIII, комната 3.2 sales@white-product.com	
Для профессионального применения Состав набора: 1. Тестовая кассета – 40 шт; 2. Флакон с буфером – 1 шт.; 3. Одноразовые капиллярные пипетки – 40 шт; 4. Набор пленочных этикеток – 2 шт; 5. Инструкция к медицинскому изделию – 1 шт.		
		
 РУ	Дата регистрации	
	  40  	

