




 LEPU MEDICAL

**SARS-CoV-2**  
**Antigen Rapid Test Kit**  
(Colloidal Gold Immunochromatography)

**SARS-CoV-2** CE  
**Antigen Rapid Test Kit**  
(Colloidal Gold Immunochromatography)

25 Tests

**Product Name: SARS-CoV-2 Antigen Rapid Test Kit  
(Colloidal Gold Immunochromatography)**



**BEIJING LEPU MEDICAL TECHNOLOGY CO., LTD.**

Building 7-1 No.37 Chaoqian Road, Changping District,  
Beijing 102200, P.R.China  
Tel: +86-10-80123964  
Email: lepuservice@lepumedical.com  
Web: en.lepumedical.com



**Lepu Medical (Europe) Cooperatief U.A.**

Abe Lenstra Boulevard 36,8448 JB, Heerenveen, The Netherlands  
Tel: +31-515-573399  
Fax: +31-515-760020

Набор реагентов «Экспресс-тест иммунохроматографический»  
для качественного обнаружения антигена коронавируса SARS-CoV-2  
в клиническом образце для диагностики in vitro (SARS-CoV-2 Antigen Rapid Test Kit  
(Colloidal Gold Immunochromatography)), серия 20CG2710X, 20CG2710X, 20CG2710X

Производитель: Бейджин Лепу Медикал  
Технологии Ко., Лтд.  
Корпус 7-1 №37 Чжаоцзян Роуд, район Чанпин,  
Пекин, Китай, 102200

Уполномоченный представитель производителя в РФ:  
ООО «РАИД ГЛОБАЛ»,  
670000, РФ, Улан-Удэ, ул. Саянских д. 13А, офис 5. Тел: +7(495) 105-99-16


Состав набора:  
1. Тест-касета (в индивидуальной упаковке с осушителем) - 25 шт.  
2. Инструкция к медицинскому изделию - 1 шт.  
3. Буферный раствор для разведения образцов - 2 шт.


Для профессионального применения


LOT 20CG2710X 22/10/2021

CE IVD


+4°C / +30°C








Use swab to collect sample




Put off adhesive



Insert swab into test tube



a Add 8 drops of buffer to well A  
b Rotate the strip, two rounds each direction



Fold left side over, and wait for 15 minutes



Reading Area



A

B



LEPU MEDICAL




(+)



Positive

(-)



Negative

# SARS-CoV-2 Antigen Rapid Test Kit

(Colloidal Gold Immunochromatography)

CE





Прошито, пронумеровано

Количество листов 6

Подпись /Жимбиев А.Ц/  
ФИО