

УТВЕРЖДАЮ
Генеральный директор
ООО «Миндрей Медикал Рус»
 **Чэнь Син**

30 декабря 2020 г.

ПО ДОВЕРЕННОСТИ ОТ 20.08.20Г
СЫСОВЕВА Н.В.



Фотографические изображения изделия:

**Набор калибраторов для количественного определения секреторного белка 4
эпидидимиса человека (HE4 Calibrators) иммунохемилюминесцентным
методом в клиническом образце на анализаторах серии CL для диагностики in
vitro**

Shenzhen Mindray Bio-Medical Electronics Co., Ltd. («Шэньчжэнь Майндрэй Био-
Медикал Электроникс Ко., Лтд.»), Китай

mindray

P/N:105-010138-00

HE4 Calibrators

Contents: C0: 1×1.2 mL, C1: 1×1.0 mL, C2: 1×1.0 mL

Калибраторы HE4
Calibradores de HE4
Calibradores de HE4
Calibratori per HE4
HE4 Kalibratörleri
Calibrateurs de HE4




SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.
Mindray Building, Keji 12th Road South, High-Tech Industrial Park,
Nanshan, Shenzhen, 518057, P. R. China

EC REP

Shanghai International Holding Corp. GmbH(Europe)
Eiffestraße 80, 20537 Hamburg, Germany

IVD



PA1202007020013401

REF HE4211

LOT 2020070200 2022-02-13

047-018035-00(1.0)

mindray
HE4 Calibrators

Product Information


Catalog No.	Package Size
HE4211	C0:1x1.2 mL/vial C1:1x1.0 mL/vial C2:1x1.0 mL/vial

Intended Use
Mindray HE4 Calibrators (HE4 CAL) are used to calibrate the Human Epididymal protein 4 (HE4) assay on Mindray CL-series Chemiluminescence Immunoassay Analyzer.

Summary
Mindray CL measurement system is composed of Mindray CL-series Chemiluminescence Immunoassay Analyzer, Mindray reagent kits, Calibrators, and controls. HE4 Calibrators contain Human Epididymal protein 4 (HE4) analyte. The mathematical relationship between the measured light emission and the known analyte concentration establishes the calibration curve. It is used to convert Relative Light Units (RLUs) measured from samples to quantitative analyte concentration.

mindray
HE4
Calibration Card

Lot# 2020070200
2022-02-13



Всего прошито, пронумеровано и скреплено
печатью и подписью: 3 листа (-ов)

Генеральный Директор Сысоева / Чэнь Син
М.П.



ПО ДОВЕРЕННОСТИ ОТ 20.08.20Г
СЫСОЕВА Н.В.