

INNO-LIA™ HIV I/II SCORE

Line immunoassay for the confirmation and discrimination of antibodies to HIV-1, HIV-1 group O, and HIV-2.



**Even a winner
can improve**

CE

Short procedure
COlor-coded reagents
REady-to-use reagents




INNOGENETICS
BIOTECHNOLOGY FOR HEALTHCARE

INNO-LIA™ HIV I/II Score


All HIV strains detected by only one strip

Approved performance



- High sensitivity for HIV-1: 100% (273/273) and HIV-2: 100% (120/120)
- High specificity on screened-negative blood donors: 96.7% (290/300)

Clear-cut results



Results based on seroconversion samples from 12 BBI panels

The clearest answer in HIV confirmation.


Software and automation



Auto-LIA™ 48 offers fully automated processing for the complete range of INNO-LIA™ tests, from sample incubation to color development.



Automated data entry using a flatbed scanner.



LiRAS™ provides scanner calibration to ensure accurate reading of the results.



LiRAS™ enables customized reports using the full "standard" report or the brief "summary" report.

Presentation

INNO-LIA™ HIV I/II Score

20 tests/kit

Art. no.

80540



INNOGENETICS
BIOTECHNOLOGY FOR HEALTHCARE

Commercial office: INNOGENETICS NV
Technologiepark 6 - B-9052 ZWIJNAARDE - BELGIUM
T +32 (0)9 290 17 70 F +32 (0)9 290 17 75

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Your Agent:

For use outside the US.


«ИННО-ЛИА ВИЧ I/II подтверждающий»

Набор реагентов для выявления антител к индивидуальным белкам ВИЧ 1 типа, ВИЧ 1 типа группы О, ВИЧ 2 типа в сыворотке или плазме крови человека методом иммунного блоттинга (линейный иммуноблот). Превосходит по чувствительности Вестерн-блот, дает в 4 раза меньше неопределенных результатов.



Производство INNOGENETICS N.V., Бельгия.

- Полный набор рекомбинантных белков и синтетических пептидов ВИЧ-1, ВИЧ-1 группы О и ВИЧ-2 (sgp120, gp41, p31, p24, p17, sgp105, gp36) позволяет охарактеризовать антитела у пациента против белков ВИЧ различных субтипов;
- Все антигены фиксированы на полоске на четко определенных позициях;
- Возможно дифференциальное подтверждение на одной полоске ВИЧ-1 и ВИЧ-2 сероконверсии;
- Наличие внутренних контрольных полос;
- Готовые к использованию жидкие реагенты.

Простая процедура анализа:



Структура тест-полоски:



Набор рассчитан на 20 определений, включая контрольные.
Срок годности набора 12 месяцев.

Handwritten signature

