

Фотографии медицинского изделия:

«Набор инструментов стоматологических полировальных KENDA»

производства KENDA AG
Kanalstrasse 18, 9490 Vaduz, Liechtenstein/ Лихтенштейн

Генеральный директор ООО «МедТестКонсалтинг»

(должность)

ЭЛЕКТРОННАЯ КОПИЯ ОФОРМЛЕНА ЦЕНТРОМ ИСПЫТАНИЙ И РЕГИСТРАЦИИ МЕДИЦИНСКИХ ИЗДЕЛИЙ

(подпись и печать Заявителя)



В. К. Бунзя
(И.О.Фамилия)

KENDA
DENTAL POLISHERS
KENDA NOBILIS
1-step diamond polishing system for composite restorations and ceramics
KENDA AG, 11 9480 Vöcklabruck
www.kenda-dental.com



REF 0300.012



CE 0120

LOT 9.7143

2022-07-17

135C
275F

* EKEN03000121/SS20717229.7143I *

12 KENDA Nobilis Diamond Polishers for
Composite and Ceramic, RA

max 10000 rpm, use water spray, autoclavable



REF 0350.012

LOT 9.7144

HIBC



CE 0120

2022-06-26

135C
275F



WMOS

* + EKEN03500121/\$\$20821229.7144K*

12 KENDA Unicus Diamond Polishers for

Composite and Ceramic, RA

max 10000 rpm, use water spray, autoclavable





Manufactured by

KENDA AG
INDUSTRIEZONE NEUGUT
KANALSTRASSE 18
LI-9490 VADUZ
LIECHTENSTEIN

www.kenda-dental.com

REF 0500.012

HIBC



CE 0120

LOT 9.7145

2022-07-17

135C
275F

* KEND05000121/5520717229.7145M*

12 KENDA Complete Polishing System for
Composite/Compomer + Gold/Amalgam, RA

max 1000 rpm, use water spray, autoclavable



REF 0380.012


LOT 9.7146



CE 0120



2022-07-17



135C
275F

1 KEND03800121/SS20717229.7146T

12 KENDA Maximus Diamond Polishers for
Composite and Hybrid-Ceramic, RA

max 10000 rpm, use water spray, autoclavable

REF 0400.012

LOT 9.7147

HIBC



CE 0120



2022-04-25



135C
275F



HM05

* EKFN04000121/\$\$20821229.7147J*

**12 KENDA Deluxe Diamond Polishers for
Composite and Hybrid-Ceramic, RA**

max 10000 rpm use water spray, autoclavable




Фото 11. VI. Набор KENDA HYBRID REF 790.012 (полир чашевидный малый светло-серый - 1 шт., полир чашевидный средний светло-серый - 1 шт., полир
 000 ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU 1 шт., полир копьевидный большой светло-серый - 1 шт., полир дисковидный светло-серый - 1 шт., полир чашевидный
 большой светло-серый - 1 шт., полир чашевидный малый бледно-розовый - 1 шт., полир чашевидный средний бледно-розовый - 1 шт., полир копьевидный малый

REF 790.012



MIBC

CE 0120

LOT 9.7148

2022-07-17


135C
275F

...KEN7900121/4520717229.7148.*

12 KENDA Polishers for

Hybrid Composites, RA

max 10000 rpm, use water spray, autoclavable



REF 791.006 HP



FKEN7910061/59.71490

6 KENDA Polishers for
Hybrid Composites, HP

max 10000 rpm, max. contact pressure 2N

LOT 9.7149

2022-08-21



Manufactured by

KENDA AG
INDUSTRIEZONE NEUCHÂT
KAMALSTRASSE 18
LI-9490 VADUZ
LIECHTENSTEIN

www.kenda-dental.com











REF 933.F.006

LOT 11933

HIBC



CE 0120

2022-07-17


135C
275F

* EKFN933F0061/552080922119336*

6 KENDA C.G.I. Polishers f. Composites

Compomers and Glass-ionomers, RA

max 10000 rpm, use water spray, autoclavable



REF 0266.006

LOT 12003



CE 0120



2022-07-17



135C
275F

REF 0266.006 / 55207172212003K*

5 KENDA Zircovis Diamond Hi-Shine Pol.

optimised for Zirconia, RA

max 10'000 rpm. use water spray, autoclavable



REF 0250.012

WIBC



CE 0120

LOT 9.7152

2022-07-17

135C
275F

* KEN02500121/\$\$20717229.7152M*

12 KENDA Zircovis Diamond Polishers

optimised for Zirconia, RA

max. 10000 rpm, use water spray, autoclavable



REF 0200.012

HIBC




CE 0120

LOT 9.7153



2021-03-04



135C
275F



ИМ05

* * EKFN02000121/\$\$20821229.7153E*

12 KENDA All Ceramic Diamond Polishers

for all ceramic restorations, RA

max. 10000 rpm, use water spray, autoclavable



REF 780.012



CE 0120

LOT 9.7154



2022-07-17



135C
275F

* KI N7800121 / S320717229.7154W*

12 KENDA Diagum Polishers

for Ceramic, RA

max. 10000 rpm, use water spray, autoclavable



REF 781.006 HP



5 KENDA Diagam Polishers
for Ceramic, HP

Lot 9.7155

2021-06-30





REF 5100.012



CE 0120

LOT 9.7156

2022-07-17

135C
275F

FKEN51000121/\$\$20717229.7156P

12 KENDA Propy Polishers
for Prophylaxis, RA

max. 10000 rpm use water spray. autoclavable



REF 0600.012



CE 0120

LOT 9.7157



2022-07-17

135C
275F

* (REF 0600.012 / 5520717229.7157Q)*

12 KENDA Planus Acrylic Denture
Polishers, RA

тип: 12 штук в коробе, автоклавируемые





Caution:

Max. speed 15'000 RPM

Max. contact pressure 2N

REF 6700.006



* EKEMA7000061/8870821229.7158*

6 KENDA Planus Polishers for Acrylic

Dentures, coarse, HP

max. 15'000 rpm, acceleration, max. contact pressure 2N

LOT 9.7158

2022-05-16

135C
275F



10M05



REF 700.012

INIBL



CE 0120

LOT 9.7159

2022-03-29

135C
275F



HM05

* FKEN7000121/\$\$20821229.7159P*

12 KENDA Polishers

for Amalgam, RA

max. 10'000 rpm, use water spray, autoclavable



REF 740.012

LOT 9.7160



CE 0120



2022-06-22



135C
275F





HM05

* EKEN7400121/\$\$20821229.7160L*

12 KENDA Polishers for metals and
titanium alloys, RA

max 10'000 rpm, use water spray, autoclavable





REF 3006.005FG LOT 11940 Rx only



max. 20'000 rpm

Min 36°F/2°C
Max 82°F/28°C



2022-07-17



Caution consult
instructions before use



See instructions
for use




0120

* + EKEN3006005FG1/\$\$2071722119408*

5 KENDA Hi-Shine Polishers
for Amalgam, FG



Manufactured by KENDA AG • Kanalstrasse 18
LI-9490 Vaduz • Liechtenstein • www.kenda-dental.com



REF Z-000.006




* EKFNZD000061/SS20821229.71627*

6 KENDA Z-Series Diamond Polishers
for all types of Ceramic, HP
max. 10 min. rpm, max. contact pressure 2N

LOT 9.7162

2021-07-21






REF 626.004

LOT 9.7163

 2022-08-21

HIBC



* FKFN6260041/5820821229.7163S*



ИМ05

4 KENDA HighEnd Diamond Power Discs

for all types of Ceramics, HP

max 10000 rpm, max contact pressure 2N





A



Б



В

 **Caution:**

Max. speed 15'000 RPM

Max. contact pressure 2N

REF 9008.006



LOT 11818

2022-08-21



* EKEN90080061/55208212211818W*
6 KENDA Polishers for Acrylic
Dentures, coarse, HP
max 15'000 rpm, max. contact pressure 2N



HE-V104.5


Caution: Read instructions
www.kenda-dental.com

Wear eye protection


Wear breathing protection

Storage temperature
Max 2°C / 34°F
Min 28°C / 82°F

Keep away from
light







 **Caution:**

Max. speed 15'000 RPM

Max. contact pressure 2N

REF 8508.006



LOT 11852

2022-08-21



6 KENDA Polishers for Acrylic
Dentures, extra-fine, HP
max. 15'000 rpm, max. contact pressure 2N



 **Caution:**

Max. speed 10'000 RPM

Max. contact pressure 2N

REF Y-000




LOT 9.7164


 2022-08-21



* KENDAY0001/6620821229.71642*

3 KENDA Yellow Line Polishers for
Acrylic Dentures, medium, HP
max. 10'000 rpm, max. contact pressure 2N





KENDA
DENTAL POLISHERS
ZIRCO 1 by KENDA
diamond polishing



KENDA
DENTAL POLISHERS
ZIRCO 1 by KENDA

REF 0230.006

LOT 9.7166



CE 0120

2022-07-17

135C
275F

* KEN02300061/SS20717229.7166S*

6 ZIRCO 1 by KENDA 1-step Diamond Pol.

for zirconia and all ceramics, RA

max 15000 rpm, use water spray, autoclavable

Прошито, пронумеровано и
скреплено печатью

58 листа(ов)

Генеральный директор
ООО «МедТестКонсалтинг»

В.К. Бунзя

