


«УТВЕРЖДАЮ»
Генеральный директор
ООО «Миндрей Медикал Рус»


Сюн Вэй
по доверенности от 31.05.21 г.
СЫСЬЕВА Н.В.
«19» октября 2021 г.

ФОТОГРАФИЧЕСКИЕ ИЗОБРАЖЕНИЯ МЕДИЦИНСКОГО ИЗДЕЛИЯ

Реагент лизирующий LC LYSE для анализатора автоматического
гематологического BC-5390 CRP

Реагент лизирующий LC LYSE для анализатора автоматического гематологического BC-5390 CRP

Hematology Analyzer Reagent

LC LYSE

The LC Lyse is intended to be used with Mindray BC-5390 CRP/BC-5180 CRP/BC-5310 CRP/BC-5100 CRP/BC-7500[B] CRP/BC-7500[N] CRP/BC-7500[R] CRP/BC-7500[NR] CRP/BC-7300[B] CRP/BC-7300[N] CRP/BC-7300[R] CRP/BC-7300[NR] CRP/BC-7500[B] CS/BC-7500[N] CS/BC-7500[R] CS/BC-7500[NR] CS Auto Hematology Analyzer and CRP-M100 Specific Protein Analyzer for dissolving RBCs before sample analysis.

Active Ingredients :

Morpholinopropanesulfonic acid buffer..... 0.60%.


Storage and Stability :

The LC lyse can be stored for 12 months at 2°C ~ 8°C.

The LC lyse should be used at the ambient temperature of 15°C ~ 32°C; once opened, the lyse may remain valid for 30 days.



0ED202101170100191



LOT 2021011701 2022-01-16 mindray

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.
Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan,
Shenzhen, 518057, P.R. China
Tel: +86 755 81888998 Fax: +86 755 26582680 Website: <http://www.mindray.com>
Shanghai International Holding Corp. GmbH(Europe)
Eiffelstraße 80, 20537 Hamburg, Germany 047-011805-00(8.0)


mindray

LC LYSE

200 mL

0093

LOT 2021011701  **2022-01-16**



00013979165320203758

IVD



SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.

Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan,
Shenzhen, 518057, P.R. China


Tel: +86 755 81888998

Fax: +86 755 26582680

Website: <http://www.mindray.com>

EC REP

Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, Germany



mindray

047-011819-00(5.0)

The LC Lyse is intended to be used with Mindray BC-5390 CRP/BC-5180 CRP/BC-5310 CRP/BC-5100 CRP/BC-7500[B] CRP/BC-7500[N] CRP/BC-7500[R] CRP/BC-7500[NR] CRP/BC-7300[B] CRP/BC-7300[N] CRP/BC-7300[R] CRP/BC-7300[NR] CRP/BC-7500[B] CS/BC-7500[N] CS/BC-7500[R] CS/BC-7500[NR] CS Auto Hematology Analyzer and CRP-M100 Specific Protein Analyzer for dissolving RBCs before sample analysis.

Active Ingredients :

Morpholinopropanesulfonic acid buffer..... 0.60%

Storage and Stability :

The LC Lyse can be stored for 12 months at 2°C ~ 8°C.

The LC Lyse should be used at the ambient temperature of 15°C ~ 32°C; once opened, the Lyse may remain valid for 30 days.

Warnings :

Do not inhale and or ingest. In case of inhalation or ingestion, seek medical advice immediately.

Avoid skin contact. In case of skin contact, flush the affected area with plenty of water immediately.

Avoid eye contact. In case of eye contact, flush the affected area with plenty of water and seek medical advice immediately.

Precautions :

Do not use product past expiration date.

Do not use the product when it is frozen(either partially or completely).

Dispose of waste product, residual product and contaminated packaging in accordance with government regulations.

Макет русскоязычной маркировки упаковки

Реагент лизирующий LC LYSE для анализатора автоматического гематологического ВС-5390 CRP

Производитель: Шэньчжэнь Майндрэй Био-Медикал Электроникс Ко., Лтд.
(Shenzhen Mindray Bio-Medical Electronics Co., Ltd.)
Mindray Building, Keji 12th Road South, High-Tech Industrial Park,
Nanshan, 518057 Shenzhen, People's Republic of China

Уполномоченный представитель в РФ: ООО «Миндрей Медикал Рус», Россия, 129090, г. Москва,
Олимпийский проспект, д.16, стр.5, антресоль 4, помещение I,
ком.7, 11А

Тел. +7 (499) 553-60-36, E-mail: info.ru@mindray.com

Регистрационное удостоверение № от



схематическое изображение
6 листов

