

Фотографические изображения общего вида на медицинское изделие

Калибраторы для обеспечения правильности количественного определения N-терминального фрагмента натрийуретического пропептида В-типа (NT-proBNP) иммунохемилюминесцентным методом в сыворотке и плазме крови человека на иммунохимических анализаторах ARCHITECT «NT-proBNP Калибраторы (Alere NT-proBNP for ARCHITECT Calibrators)»

(размером не менее 18 x 24 сантиметра)




Alere NT-proBNP for ARCHITECT Calibrators

		pg/mL	pmol/L
CAL A	1 x 4.0 mL	0.0	0.0
CAL B	1 x 4.0 mL	120.0	14.2
CAL C	1 x 4.0 mL	220.0	26.0
CAL D	1 x 4.0 mL	345.0	40.7
CAL E	1 x 4.0 mL	4100.0	483.8
CAL F	1 x 4.0 mL	41,000.0	4838.0



UNTIL FIRST USE



AFTER THAW

CONTAINS: AZIDE




EC REP


DISTRIBUTED BY
Abbott Laboratories

Axis-Shield Diagnostics Limited
Luna Place, The Technology Park,
Dundee DD2 1XA
United Kingdom
+44-1382-422000

MDSS GmbH
Schiffgraben 41,
30175 Hannover, Germany
Phone: (+49)-511-5262-6630
Fax: (+49)-511-5262-6633

PRODUCT OF GB







(01) 05055845400510 (17) 210912
(10) 902936107 (240) 2R1002
2021-09-12
902936107
LOT
R01
www.alere.com
Alere NT-proBNP Calis
2R10-02

Alere NT-proBNP for ARCHITECT Calibrators

		pg/mL	pmol/L
CAL A	1 x 4.0 mL	0.0	0.0
CAL B	1 x 4.0 mL	120.0	14.2
CAL C	1 x 4.0 mL	220.0	26.0
CAL D	1 x 4.0 mL	345.0	40.7
CAL E	1 x 4.0 mL	4100.0	483.8
CAL F	1 x 4.0 mL	41,000.0	4838.0

20°C
UNTIL FIRST USE
AFTER THAW

CONTAINS AZIDE
!
DO NOT
DISTRIBUTED BY
Abbott

Abbott
Abbott Laboratories
Abbott Park, 100 Abbott Park, IL 60014-1000
USA
Tel: (708) 462-1000
Fax: (708) 462-1001
Email: info@abbott.com
Website: www.abbott.com

ARCHITECT
Alere NT-proBNP for ARCHITECT Calibrators
Created March 2020

en
Alere NT-proBNP for ARCHITECT
2R10
ABBL634/R01
S2R120

INTENDED USE

The Alere NT-proBNP for ARCHITECT Calibrators are for the calibration of the ARCHITECT iSystem when used for the quantitative determination of N-terminal pro B-type natriuretic peptide (NT-proBNP) in human serum and plasma.

Refer to the reagent package insert and the ARCHITECT System Operations Manual for additional information.

CONTENTS

6 Bottles (4.0 mL each) of Alere NT-proBNP for ARCHITECT Calibrators. Calibrators A-F contain Tris buffer with protein (bovine) stabilizers. Calibrators B-F also contain different concentrations of recombinant NT-proBNP. Preservatives: ProClin 300 and sodium azide.

The calibrators are at the following concentrations:

Response	
P302+P362 P333+P313	IF ON SKIN: Wash with plenty of water. If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

* Not applicable where regulation EC 1272/2008 (CLP) has been implemented.

Follow local chemical disposal regulations based on your location along with recommendations and content in the Safety Data Sheet to determine the safe disposal of this product.

For the most current hazard information, see the product Safety Data

Всего прошито, пронумеровано
и скреплено печатью

4 листов



Подпись

Макет маркировки на медицинское изделие

Калибраторы для обеспечения правильности количественного определения N-терминального фрагмента натрийуретического пропептида В-типа (NT-proBNP) иммунохемилюминесцентным методом в сыворотке и плазме крови человека на иммунохимических анализаторах ARCHITECT «NT-proBNP Калибраторы (Alere NT-proBNP for ARCHITECT Calibrators)»



REF 2R10-02

IVD



Калибраторы для обеспечения правильности количественного определения N-терминального фрагмента натрийуретического пропептида В-типа (NT-proBNP) иммунохемилюминесцентным методом в сыворотке и плазме крови человека на иммунохимических анализаторах ARCHITECT «NT-proBNP Калибраторы (Alere NT-proBNP for ARCHITECT Calibrators)»

В составе:

- Калибратор А (Alere NT-proBNP Calibrator A) – 1 x 4,0 мл
- Калибратор В (Alere NT-proBNP Calibrator B) – 1 x 4,0 мл
- Калибратор С (Alere NT-proBNP Calibrator C) – 1 x 4,0 мл
- Калибратор D (Alere NT-proBNP Calibrator D) – 1 x 4,0 мл
- Калибратор E (Alere NT-proBNP Calibrator E) – 1 x 4,0 мл
- Калибратор F (Alere NT-proBNP Calibrator F) – 1 x 4,0 мл
- Инструкция по применению

Производитель: Аксис-Шилд Диагностикас Лимитед, Axis-Shield Diagnostics Limited, Luna Place, The Technology Park, Dundee, DD2 1XA, United Kingdom (Соединенное Королевство)

Место производства: Axis-Shield Diagnostics Limited, Luna Place, The Technology Park, Dundee, DD2 1XA, United Kingdom (Соединенное Королевство)

Страна происхождения: Великобритания / РУ № ____/____ от ____/____/ Инструкция по применению на русском языке: www.corelaboratory.abbott/IFU

Серию, срок годности, условия хранения, страну происхождения и графические символы см. на оригинальной упаковке. Расшифровку графических символов см. в инструкции.

Авторизованный представитель изготовителя на территории России и других стран СНГ: Общество с ограниченной ответственностью «Эбботт Лэбораториз», 125171, Россия, г. Москва, Ленинградское шоссе, 16А, строение 1,
Тел: +7 (495) 258 42 70 (80); Факс: +7 (495) 258 42 71 (81)

ABVL635P03


www.corelaboratory.abbott/IFU

Alere NT-proBNP Cals 2R10-02 REF




Alere NT-proBNP for ARCHITECT Calibrators

		pg/mL	pmol/L
CAL A	1 x 4.0 mL	0.0	0.0
CAL B	1 x 4.0 mL	120.0	14.2
CAL C	1 x 4.0 mL	220.0	26.0
CAL D	1 x 4.0 mL	345.0	40.7
CAL E	1 x 4.0 mL	4100.0	483.8
CAL F	1 x 4.0 mL	41 000.0	4838.0



UNTIL FIRST USE



AFTER THAW

CONTAINS: AZIDE



EC REP

DISTRIBUTED BY
Abbott Laboratories

PRODUCT OF GB

Axis-Shield Diagnostics Limited
Luna Place, The Technology Park,
Dundee DD2 1XA
United Kingdom
+44-1382-422000

Abbott Rapid Dx International Limited
Parkmore East Business Park,
Ballyone,
Co. Galway, H91 VK7E,
Ireland
Tel.: +353 91 429 900


Abbott

Alere NT-proBNP
Cals

Калибратор А (Alere NT-proBNP Calibrator A):

Alere NT-proBNP for ARCHITECT

CAL A



-20°C

CONTAINS: AZIDE



IVD

REF

2R10A

4.0 mL

0.0 pg/mL

ABBL636R02

Axis-Shield Diagnostics Ltd




Exp.

LOT

Калибратор В (Alere NT-proBNP Calibrator B):

Alere NT-proBNP for ARCHITECT

CAL B



-20°C

CONTAINS: AZIDE



IVD

REF

2R10B

4.0 mL

120.0 pg/mL

ABBL637R02

Axis-Shield Diagnostics Ltd



Exp.

LOT

Калибратор С (Alere NT-proBNP Calibrator C):

Alere NT-proBNP for ARCHITECT

CAL C



-20°C

CONTAINS: AZIDE



IVD

REF

2R10C

4.0 mL

220.0 pg/mL

ABBL638R02

Axis-Shield Diagnostics Ltd



Exp.

LOT

Калибратор D (Alere NT-proBNP Calibrator D):

Alere NT-proBNP for ARCHITECT

CAL D



-20°C

ABBL639R02

Axis-Shield Diagnostics Ltd

CONTAINS: AZIDE



Exp.

LOT

IVD

REF 2R10D

4.0 mL

345.0 pg/mL



Калибратор E (Alere NT-proBNP Calibrator E):

Alere NT-proBNP for ARCHITECT

CAL E



-20°C

ABBL640R02

Axis-Shield Diagnostics Ltd

CONTAINS: AZIDE



Exp.

LOT

IVD

REF 2R10E

4.0 mL


4100.0 pg/mL



Калибратор F (Alere NT-proBNP Calibrator F):

Alere NT-proBNP for ARCHITECT

CAL F



-20°C

ABBL641R02

Axis-Shield Diagnostics Ltd

CONTAINS: AZIDE



Exp.

LOT

IVD

REF 2R10F

4.0 mL

41,000.0 pg/mL



Всего прошито, пронумеровано
и скреплено печатью

4 листов

