



ООО «Бектон Дикинсон
Восток»
127051, Российская Федерация,
Москва, ул. Садовая-
Самотечная,
24/27, 2 этаж, помещение I,
комната 18
тел: + 7 495 775 8582
факс: + 7 495 775 8583
bd.com/ru

ПОДТВЕРЖДАЮ

Специалист по регистрации
по России и странам СНГ
ООО «Бектон Дикинсон Восток»
Борзых М. В.
«29» июня 2022 г.



**Фотографические изображения медицинского изделия для диагностики
in vitro**

Набор для подготовки цитологических препаратов для прибора для
автоматизации пробоподготовки перед проведением цитологических
исследований BD Totalys MultiProcessor



Рисунок 1. Реагент плотности, 480 мл




Рисунок 2. Реагент плотности, 1,9 л

Sigma P.O. Box 14508
St. Louis, MO 63178 USA

T-6664

Tris Buffered Saline, pH 8.0

1 Packet Lot SLCF8558

Dry powder in foil pouch; when dissolved in 1 L. deionized water, will yield 0.05M Tris, 0.138M NaCl, 0.0027M KCl, pH8.0 at 25° C. Certain pH electrodes do not give accurate results with Tris buffers.

Store at room temperature.

IRRITANT:

Irritating to eyes, respiratory system and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing.

For R&D use only. Not for drug, household, or other use.
SDS Available

Sigma P.O. Box 14508
St. Louis, MO 63178 USA

T-6664

Tris Buffered Saline, pH 8.0

1 Packet Lot SLCF8558

Dry powder in foil pouch; when dissolved in 1 L. deionized water, will yield 0.05M Tris, 0.138M NaCl, 0.0027M KCl, pH8.0 at 25° C. Certain pH electrodes do not give accurate results with Tris buffers.

Store at room temperature.

IRRITANT:

Irritating to eyes, respiratory system and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing.

For R&D use only. Not for drug, household, or other use.
SDS Available



Рисунок 5. Шприц-пипетки, 480 шт.



Рисунок 6. Центрифужные пробирки, 480 шт.



Синто, пронумеровано 6
Борис М.В.

