

Over the wire (OTW) PTA balloon catheter

# JOKER® 014

PTA balloon technology specifically designed for  
intrapopliteal interventions

Q2 Eurocor Size 2.50 x 40 mm  
JOKER® 014



PTA balloon technology specifically designed for  
intrapopliteal interventions

Q2 Eurocor Size 2.50 x 40 mm  
JOKER® 014



PTA balloon technology specifically designed for  
intrapopliteal interventions

Q2 Eurocor Size 2.50 x 40 mm  
JOKER® 014



PTA balloon technology specifically designed for  
intrapopliteal interventions

Q2 Eurocor  
JOKER® 014



CE 0434

MD 214-2540



NN 2011-67

LOT 123609

EXP 2013-07

Over the wire (OTW\*) PTA Balloon Catheter  
Over the wire (OTW\*) PTA balloon catheter  
Cathéter à ballonnet à fil guidage (OTW\*) pour PTA  
Catéter con balón PTA con guía (OTW\*)  
Catetere a pallonnetto da PTA over the wire (OTW\*)

0.014"  
(0.36 mm)



| Inflation Pressure<br>atm / bar / kPa | Balloon Diameter (mm) |      |      |      |      |
|---------------------------------------|-----------------------|------|------|------|------|
|                                       | 2.0                   | 2.5  | 3.0  | 3.5  | 4.0  |
| 4                                     | 1.88                  | 2.37 | 2.79 | 3.27 | 3.87 |
| 6                                     | 2.00                  | 2.50 | 3.00 | 3.50 | 4.00 |
| 8                                     | 2.08                  | 2.65 | 3.12 | 3.72 | 4.12 |
| 10                                    | 2.16                  | 2.70 | 3.28 | 3.87 | 4.19 |
| 12                                    | 2.25                  | 2.75 | 3.39 | 4.00 | 4.27 |
| 14                                    | 2.27                  | 2.83 | 3.49 | 4.07 | 4.35 |
| 16                                    | 2.33                  | 2.89 | 3.54 | 4.12 | 4.43 |
| 18                                    | 2.43                  | 2.94 | 3.62 | 4.26 | 4.50 |
| 20                                    | 2.50                  | 3.00 |      |      |      |


OTW 0.014 in (0.36 mm)  
Recommended Guidewire

BC 0.46 in (11.68 mm) / 3.33 mm  
Recommended Sheath compatibility

Manufacturer: Eurocor Medical GmbH, D-82027 Bad Reichenhagen, Germany. Tel: +49 (0) 89 30 00 00. Fax: +49 (0) 89 30 00 01. Email: info@eurocor-med.com

Bitte lesen Sie die Anweisung an dem  
Packaging sehr aufmerksam zu.  
18.7.11 Eurocor GmbH  
in den Dänen 6a  
D-82117 Bogen  
Tel: +49 (0) 89 30 00 00  
Fax: +49 (0) 89 30 00 01  
Email: info@eurocor-med.com





*A. B.*

Сидоров А.В.

Over the wire (OTW) PTA balloon catheter

# JOKER® 035

PTA balloon technology specifically designed for  
intraoperative interventions

Eurocor Size: 6.0x150 mm  
JOKER® 035 Lot: No. 127095



Eurocor Size: 6.0x150 mm  
JOKER® 035 Lot: No. 127095



Eurocor Size: 6.0x150 mm  
JOKER® 035 Lot: No. 127095



Eurocor  
JOKER® 035



Over the wire (OTW\*) PTA Balloon catheter  
Over the wire (OTW\*) PTA balloon catheter  
Cathéter à ballonnet à l'hydrodilatation (OTW\*) pour PTA  
Catéter con balón PTA con guía (OTW\*)  
Catetere a palloncino da PTA over the wire (OTW\*)

CE 0434


REF: 435-60150 S

EXP: 2012-02

LOT: 127095

EXP: 2014-02

0.035"  
(0.89 mm)




GW: ≤ 0.035" in (0.89 mm)  
Recommended Guidewire

BC: ≥ 6F (0.067 in / 1.67 mm)  
Recommended Sheath compatibility

| Inflation Pressure<br>PSI / bar / MPa | Balloon Diameter (mm) |      |      |      |      |
|---------------------------------------|-----------------------|------|------|------|------|
|                                       | 4.0                   | 5.0  | 6.0  | 7.0  | 8.0  |
| 4                                     | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |
| 6                                     | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |
| 8                                     | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |
| 10                                    | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |
| 12                                    | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |
| 14                                    | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |
| 16                                    | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |
| 18                                    | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |
| 20                                    | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |
| 22                                    | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |
| 24                                    | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |
| 26                                    | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |
| 28                                    | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |
| 30                                    | 4.00                  | 5.00 | 6.00 | 7.00 | 8.00 |

Supplier details: Eurocor Medical, Inc., 10000 N. 10th Ave., Suite 100, Denver, CO 80231, USA  
Phone: +1 303 751 1000 Fax: +1 303 751 1001 Email: info@eurocor.com





*Handwritten signature*

Пронумеровано, прошито и  
скреплено печатью на 2  
24 листах  
Ген. директор ООО «ФОРЛАЙФ»  
(А.В.Сидоров)  
«24» июня 2014 г.

