


Фотографические изображения
медицинского изделия «Наборы
инфузионные с/без инъекционным узлом»

Бер интрузионный Winged Infusion Set, без интрузион-
но узла (Узна Губера 19Ga x 2cm)



Без интрузионной иглы (игла Губера 22 Ga x 2.5 cm)

Winged Infusion Set		REF 0604220	
		BARD	
Gauge	22 Ga.		
Needle Length	1.00 in (2.5 cm)	Bore to Tip Length	
Length	8.0 in (20.3 cm)		
Priming Capacity	0.3 cc		

PVC TUBING DOES NOT CONTAIN DEHP

Made in Mexico by
Bard Access Systems - Salt Lake City, UT

Caution: Non-Pyrogenic unless the bag is damaged or opened.
For use, this product may be a sterile disposable. Handle with care and in accordance with standard medical practice and applicable local, state and federal laws and regulations.

Attention: Use Instructions for Use. Single Use. Not Reusable.

Lot: RERG0642 Use By: 2012-08

Assembled in Mexico
079292 / 0409R / 00223



Вингед инфузионный набор
 Winged Infusion Set, с
 инъекционным узлом (19Ga x 1,3cm и на губера)



Бер инфузионный Winged Infusion Set, с инфузионным узлом (иголка 22 га x 1,3 см)


Winged Infusion Set With Injection Site		REF 0603310	
		BARD	
Gauge:	22 Ga.		
Needle Length:	0.5 in. (1.3 cm)	Bard to Top Length	
Length:	8.0 in. (20.3 cm)		
Priming Capacity:	0.3 cc		
Does Not Contain DEHP			
Bard Access Systems 405 North 5600 West Salt Lake City, UT 84116 Phone: (801) 577-5000 www.bardaccess.com			
Attention: See Instructions For Use			
STERILE EO			
Caution: Federal (U.S.A.) law restricts sale of this device to or on the order of a physician.			
LOT		Use By	
REXB0860		2018-02	
Assembled in Mexico			



Бер инфузионный узел (Winged Infusion Set, 19G x 3.8cm)




 000 ЦИПМИ  INFO@CIRMI.RU  WWW.GOSZDRAVNADZOR.RU



Набор интрузионный Winged Infusion Set, (Игла Губера 25 Ga x 38cm)

Winged Infusion Set With Injection Site		REF 0604260
BAIRD		
Gauge:	25 Ga.	
Needle Length:	1.50 in. (3.8 cm)	Send to T121615
Length:	8.0 in. (20.3 cm)	
Priming Capacity:	0.2 cc	
Does Not Contain DEHP		
Bard Access Systems 400 North 100th Street Salt Lake City, UT 84116 Tel: 801.572.1000 www.bardaccess.com		



STERILE CO. Caution: Federal (U.S.A.) law restricts sale of this device to the order of a physician.

LOT
R1XB0661

Use By
2018-02

Authorized Medical
Device 10048-1-00042



