

ПОДТВЕРЖДАЮ

Генеральный директор
АО «Деалмед»



Алексеев А.П.
«12» мая 2023 г.

**Фотографические изображения общего вида медицинского изделия
«Дефибриллятор автоматический наружный Heart Guardian»**



Фото 1 – «Дефибриллятор автоматический наружный Heart Guardian HR-501»: основной блок




Фото 2 – Электрод для дефибрилляции Р-303



Фото 3 – Батарея непerezаряжаемая BT-303



Фото 4 – Чехол для хранения портативный HR-B1



P-303 AED PAD

HEART  GUARDIAN

Product name:

External Electric Defibrillation Pads

Model: P-303

Packing unit: 1 pair of pads per 1 pouch



Usage: Delivers electric shock to the patient,
Detect Patient's ECG signal

CE
2460



RADIAN QBIO Co., Ltd.

#1609~1611, HALLA-SIGMA VALLEY,
53, Gasan digital 2-ro, Geumcheon-gu,
Seoul, Korea | www.radian.co.kr
Made in KOREA



S.B. PHARMA GMBH

Max-Planck Str. 39a D-50858, Köln
Germany

Tel) +49 (0)2234 988 1521



Disposable,
Do not re-use



⚠ Caution and Warning

1. This pouch bag contains disposable pads. Do not reuse.
2. Do not open the package before use.
3. When the package was opened, the pads should be used immediately and replace the pads every 24 hours for long term use.
4. Do not fold or apply excessive force to the package.
5. limitation machine: this model (P-303) can only be used in conjunction with Heart Guardian HR-501.
6. Defibrillation shock is very high level of electric voltage. Do not touch the pads or patients when shock is delivered.

LOT

RD210222A



2023-05-21

**2020.06.03 / REV003



Technology for Life

RADIAN QBIO

Фото 5 – Упаковка электрода для дефибрилляции P-303



Фото 6 – Маркировка дефибриллятора автоматического наружного Heart Guardian HR-501

