

Urine Sediment Analyzer Reagent



Product	Sheath Fluid
Use	Used to wrap sample, form sheath flow, and to count visual components in urine
Applicable Analyzer	FUS-100 Urine Sediment Analyzer
Package	10L / Barrel, 20L / Barrel
Product	Diluent
Use	Used to dilute urine sample
Applicable Analyzer	FUS-100 Urine Sediment Analyzer
Package	500ml/Bottle×4
Product	Detergent
Use	Used to clean the pipeline and counting chamber
Applicable Analyzer	FUS-100 Urine Sediment Analyzer
Package	500ml/Bottle×4
Product	Focus
Use	Used to determine the location of focal surface of the particle camera system
Applicable Analyzer	FUS-100 Urine Sediment Analyzer
Package	125ml/Bottle
Product	Negative Control
Use	Used for Urine Sediment Quality Control
Applicable Analyzer	FUS-100 Urine Sediment Analyzer
Package	125ml/Bottle
Product	Positive Control
Use	Used for Urine Sediment Quality Control
Applicable Analyzer	FUS-100 Urine Sediment Analyzer
Package	125ml/Bottle
Product	Standard Solution
Use	Used for Urine Sediment Calibration
Applicable Analyzer	FUS-100 Urine Sediment Analyzer
Package	125ml/Bottle