

ULTRAVERSE[®]018

PVA Balloon Catheter System

↔ Balloon Length **2 - 6 cm**

■ OP Operating Pressure ■ RBP Rated Burst Pressure

ATM	2.0mm	2.5mm	3.0mm	3.5mm	4.0mm	5.0mm	6.0mm	7.0mm	8.0mm	9.0mm
6	2.00	2.50	3.00	3.50	4.00	5.00	6.00	7.00	8.00	9.00
7	2.01	2.52	3.05	3.56	4.07	5.08	6.08	7.11	8.10	9.14
8	2.02	2.55	3.09	3.61	4.13	5.15	6.15	7.22	8.21	9.27
9	2.03	2.57	3.13	3.66	4.19	5.22	6.22	7.32	8.31	9.39
10	2.04	2.60	3.17	3.71	4.25	5.28	6.30	7.42	8.41	9.52
11	2.05	2.62	3.20	3.75	4.30	5.35	6.37	7.51	8.51	9.64
12	2.06	2.64	3.23	3.79	4.36	5.41	6.44	7.59	8.60	
13	2.07	2.67	3.26	3.83	4.41	5.47	6.51			
14	2.08	2.69	3.29	3.87	4.45	5.53	6.56			
15	2.09	2.72	3.32	3.90	4.50					
16	2.10	2.74	3.35	3.93	4.54					


BARD | PERIPHERAL VASCULAR

Manufacturer
Bard Peripheral Vascular, Inc.
1625 West 3rd Street
Tempe, AZ 85281
USA

Tel: 1-800-894-9515
1-800-321-4254
Fax: 1-480-965-7062
1-800-440-5376
www.bardpv.com


PKY232100 Rev 1 1/2008

ULTRAVERSE® 018



Contents

PTA Balloon
Dilatation Catheter



SL Shaft Length: 75 cm

REF Catalog Number U87532

LOT Lot Number 50057976

Stylet Included
Use By 2016-10

BALLOON DIAMETER 3.0 mm CATHETER LENGTH 2 cm OP DISTANCE 6 cm RBP DISTANCE 16 cm BALLOON WALL THICKNESS 3.9F CATHETER WALL THICKNESS 0.018 in / 0.5 mm BALLOON LENGTH 4F



Attention, See Instructions For Use.



Single Use.



Do Not Re-Sterilize.

Non-Pyrogenic

Non-Pyrogenic

STERILE EO

Sterilized By Using Ethylene Oxide.



Keep Dry.

Rx Only



Protect From Heat.



Do Not Use If Package Is Damaged Or Opened.



U.S. and Foreign Patents Pending.
Copyright © 2005 C. R. Bard, Inc. All Rights Reserved. Printed in Ireland.



Bard and Ultrasense are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Manufacturer:

Bard Peripheral Vascular, Inc.
1625 West 3rd Street
Tempe, AZ 85281
USA

Tel: 1-480-854-9515
1-800-321-4254

Fax: 1-480-966-7062
1-800-440-5376

www.bardpv.com
Made in Ireland

BARD | PERIPHERAL VASCULAR

CE
0000
Bard Peripheral Vascular
000000 000000 000000



ULTRAVERSE[®]

BALLOON PRESSURE
(atm)


6

ULTRAVERSE RX PTA Dilatation Catheter

BARD

PERIPHERAL VASCULAR

Contents



3.5mm
x
80mm

	BALLOON DIAMETER	3.5mm
	BALLOON LENGTH	80mm
	SL SHAFT LENGTH	150cm
	RECOMMENDED GUIDEWIRE	.014"
	RECOMMENDED INTRODUCER	4F
	RECOMMENDED GUIDE CATHETER	6F

BALLOON PRESSURE (atm)	BALLOON (mm)
6	NP 3.50
7	3.55
8	3.59
9	3.63
10	3.68
11	3.72
12	3.76
13	3.80
14	3.84
15	RBP 3.88

1 atm = 101.3 kPa

REF

Catalogue Number

U41503H8RX

LOT

Lot Number

50060008

Use By

2014-12

ULTRAVERSE RX

BARD PERIPHERAL VASCULAR

REF U41503H8RX LOT 50060008 2014-12



*41503H8RX



*50121450360008

ULTRAVERSE RX

BARD PERIPHERAL VASCULAR

REF U41503H8RX LOT 50060008 2014-12



*41503H8RX



*50121450360008

ULTRAVERSE RX

BARD PERIPHERAL VASCULAR

REF U41503H8RX LOT 50060008 2014-12



*41503H8RX




*50121450360008

ULTRAVERSE RX

BARD PERIPHERAL VASCULAR

REF U41503H8RX LOT 50060008 2014-12



*41503H8RX

1 Consult Instructions For Use.

2 Single Use.

3 Do Not Re-sterilize.

Rx only

Keep Away From Heat

Do Not Contain Latex

Keep Dry

Sterilized Using Ethylene Oxide

Do Not Use if the Product Sterilization Barrier or its Packaging is Compromised.

Non-Pyrogenic

Bard and UltraVerse are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 2012 C. R. Bard, Inc. All Rights Reserved. Printed in the USA.

ULTRAVERSE RX PTA Dilatation Catheter

Cathéter de dilatation pour PTA
PTA Dilatation Catheter
Catheter per dilatation per PTA
Catheter de dilatation pour PTA
PTA Dilatation Catheter
Catheter de dilatation pour PTA

Indicaciones: Tratamiento de PTA
PTA Dilatation Catheter
Catheter per dilatation per PTA
Catheter de dilatation pour PTA
PTA Dilatation Catheter
Catheter de dilatation pour PTA

PTA Dilatation Catheter
Catheter per dilatation per PTA
Catheter de dilatation pour PTA
PTA Dilatation Catheter
Catheter de dilatation pour PTA

BARD PERIPHERAL VASCULAR

CE

ULTRAVERSE RX
REF: U41503H8RX LOT: 50060008 2014-12
ULTRAVERSE RX
REF: U41503H8RX LOT: 50060008 2014-12

ULTRAVERSE RX

Catheter with PTA
PTA Dilatation Catheter
Catheter for dilation of PTA
Catheter for dilation of PTA
PTA Dilatation Catheter
Catheter for dilation of PTA

BARD | PERIPHERAL VASCULAR

Manufacturer:
Bard Peripheral Vascular, Inc.
1425 West 3rd Street
Tempe, AZ 85281
USA
Tel: 1-800-884-4575
Fax: 1-480-884-7060
www.bardiv.com



CE
0066
Bard Limited
Ireland
Bard Ltd
UK

PK1293300 Rev. 01/07/12

ULTRAVERSE[®] R

PTA Dilatation Catheter

With
BA-CROSS[™]
Polymer Layer
Hydrophilic
Coating

ULTRAVERSE® 014

PFA Balloon Dilatation Catheter

↔ Balloon Length 8 - 30 cm

■ OP Operating Pressure ■ RBP Rated Burst Pressure

ATM	1.5mm	2.0mm	2.5mm	3.0mm	3.5mm	4.0mm	5.0mm
6	1.50	2.00	2.50	3.00	3.50	4.00	5.00
7	1.52	2.01	2.52	3.05	3.56	4.07	5.08
8	1.53	2.02	2.55	3.09	3.61	4.13	5.15
9	1.55	2.03	2.57	3.13	3.66	4.19	5.22
10	1.56	2.04	2.60	3.17	3.71	4.25	5.28
11	1.58	2.05	2.62	3.20	3.75	4.30	5.35
12	1.60	2.06	2.64	3.23	3.79	4.36	5.41
13	1.61	2.07	2.67	3.26	3.83	4.41	5.47
14	1.63	2.08	2.69	3.29	3.87	4.45	
15	1.65	2.96	2.72	3.32	3.96	4.58	
16	1.66						


BARD | PERIPHERAL VASCULAR

Manufacturer
Bard Peripheral Vascular, Inc.
1625 West 3rd Street
Tempe, AZ 85281
USA


Tel: 1-480-894-9515
1-800-321-4254
Fax: 1-480-966-7062
1-800-445-5376
www.bardpv.com

PK1291300 Rev. D 06/12

ULTRAVERSE 014



SL Shaft Length: 150cm



Stylet Included

REF Catalog Number **U4150210**

LOT Lot Number **50052044**

Use By **2016-05**



BALLOON DIAMETER



BALLOON LENGTH

OP

OPERATING PRESSURE

RBP

RAT BALLOON PRESSURE



SHAFT DIAMETER



BORE DIAMETER



BORE LENGTH

2.0 mm

10 cm

6 Atm

15 Atm

2.8-3.6F

0.014 in / 0.4 mm

4F



Attention. See Instructions For Use.



Single Use.



Do Not Re-Sterilize.

NON PYROGENIC Non-Pyrogenic.

STERILE EO

Sterilized By Using Ethylene Oxide.



Keep Dry.

Rx Only



Protect From Heat.



Do Not Use If Package Is Damaged Or Opened.



U.S. and Foreign Patents Pending.
Copyright © 2009 C. R. Bard, Inc. All Rights Reserved. Printed in the USA.



Bard and Ultrasound are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.



Manufacturer:

Bard Peripheral Vascular, Inc.
1625 West 3rd Street
Tempe, AZ 85281
USA

Tel: 1-800-894-9515
1-800-321-4254
Fax: 1-800-966-7062
1-800-440-5376

www.bardpv.com
Made in Ireland

BARD

PERIPHERAL VASCULAR



EC REP
Bard Limited
Crawley
RH11 9BP, UK



~1801341902100~



~18051850052044~

