



Фотографические изображения медицинского изделия

«Цитофлуориметр проточный BD FACSLytic для диагностики in vitro
с принадлежностями»

Производитель: Becton, Dickinson and Company, BD Biosciences
(«Бектон, Дикинсон энд Компани, БД Байосайенсиз»), США



Специалист по регистрации
медицинских изделий
ООО «Бектон Дикинсон Восток»
Вилкова Д.И.
М.П.





«Цитофлуориметр проточный BD FACSLyric для диагностики in vitro», в вариантах исполнения, с принадлежностями

Варианты исполнения: «Цитофлуориметр проточный BD FACSLyric (2 лазера, 4 флуоресцентных параметра)»; «Цитофлуориметр проточный BD FACSLyric (2 лазера, 6 флуоресцентных параметра)»; «Цитофлуориметр проточный BD FACSLyric (3 лазера, 8 флуоресцентных параметра)»; «Цитофлуориметр проточный BD FACSLyric (3 лазера, 10 флуоресцентных параметра)»; «Цитофлуориметр проточный BD FACSLyric (3 лазера, 12 флуоресцентных параметра)»

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU















Комплект резервуаров









Крышка для резервуара






Крышка для резервуара с уплотнителем






Набор предупреждающих наклеек (1)





Набор предупреждающих наклеек (2)
















Резервуар малый






Инструмент для снятия уплотнительного кольца






Адаптер для пробирок, 15 мл



Адаптер для пробирок, 50 мл




Адаптер для пробирок Эппендорф



Office Home & Business

1 PC/Mac
2019

PRODUCT KEY INSIDE – NO DISC



Обеспечение программное MS Office




BD FACSuite™

Flow Cytometry Software

Application USB

v1.4

700023782




For Research Use Only.
Not for use in diagnostic or therapeutic procedures.

BD, the BD Logo and FACSuite are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners.
© 2019 BD. All rights reserved.

23-21572-00

Обеспечение программное BD FACSuite Software





BDFACSuite™ Clinical

700023910

BD Tritest™ 3-Color Assays

BD Tritest™ CD3/CD4/CD45

BD Tritest™ CD3/CD4/CD45 with BD Trucount™ Tubes

BD Tritest™ CD4/CD8/CD3 with BD Trucount™ Tubes


Assay USB
v1.1

BD, the BD Logo, FACSuite, Tritest and Trucount are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2019 BD. All rights reserved.


IVD CE
23-21565-00


USB 3






Модуль программный BD Multitest 4-Color Assay Module







Модуль для подачи проб универсальный BD FACS Universal Loader




Штатив на 30 пробирок (12 x 75 мм)



Штатив на 40 пробирок (12 x 75 мм)




Планшет 384-луночный




Поддон 96-луночный



Поддон U-образный 96-луночный



Капельная ловушка, левая и Капельная ловушка, правая



Пробирки тестовые





Мышь компьютерная



Диск DVD-R



BD™ Firmware Installer

REF 664238

For BD FACSLytic™ flow cytometer

v1.0.40.8



For Use by BD service personnel only

BD, the BD Logo and FACSLytic are trademarks of Becton, Dickinson and Company or its affiliates.
All other trademarks are the property of their respective owners. © 2019 BD. All rights reserved.

23-22285-00

USB-накопитель с первичной прошивкой



Монитор 23"



Монитор 27"



Считыватель штрих-кодов



Набор резервуаров для увеличения времени непрерывной работы

Макеты заводской маркировки Цитофлуориметра проточного BD FACSLytic

Model:  **BD FACSLytic™ Flow Cytometer**

REF 651164

Complies with FDA performance standards for laser products except for deviations pursuant to laser notice No. 50, dated June 24, 2007. Class 1 Laser Product per IEC/EN 60825-1:2014. This device complies with Part 15 of the FCC Rules. Operation is subjected to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This Class A digital apparatus meets all requirements of the Canadian Interference Causing Equipment Regulations.

Appareil à laser de classe 1, IEC/EN 60825-1:2014. Cet appareil numérique de la classe A respecte toutes les exigences du Règlement sur le matériel brouilleur du Canada.



2022-12-01

SN:

R651164000176

VAC-Hz:

100-240 VAC (50-60 Hz)

Total Power:

200 W



The BD FACSLytic™ flow cytometer is covered by one or more of the following US patents: 6,510,007; 7,787,197; 6,683,314; 6,809,804; 7,129,505; 6,897,954; 8,528,427



Becton, Dickinson and Company, BD Biosciences, 2350 Qume Drive, San Jose, CA 95131 USA




Becton Dickinson Ireland Ltd., Donore Road, Drogheda, Co. Louth, A92 YW26, Ireland



BD Switzerland Sàrl, Terre Bonne Park - A4, Route de Crassier 17, 1262 Eysins, Switzerland

BD, the BD Logo, and FACSLytic are trademarks of Becton, Dickinson and Company or its affiliates.
© 2021 BD. All rights reserved.

Made in USA
23-20838-02



(01)00382906511645(11)221201(21)R651164000176

Макет заводской маркировки Цитофлуориметра проточного BD FACSLytic
(пример для варианта исполнения: «Цитофлуориметр проточный BD FACSLytic
(2 лазера, 4 флуоресцентных параметра)»)

Model:  **BD FACSLytic™ Flow Cytometer**

REF 651165

Complies with FDA performance standards for laser products except for deviations pursuant to laser notice No. 50, dated June 24, 2007. Class 1 Laser Product per IEC/EN 60825-1:2014. This device complies with Part 15 of the FCC Rules. Operation is subjected to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This Class A digital apparatus meets all requirements of the Canadian Interference Causing Equipment Regulations.

Appareil à laser de classe 1, IEC/EN 60825-1:2014. Cet appareil numérique de la classe A respecte toutes les exigences du Règlement sur le matériel brouilleur du Canada.



2022-11-23

SN:

R651165000266

VAC-Hz:

100–240 VAC (50–60 Hz)

Total Power:

200 W



The BD FACSLytic™ flow cytometer is covered by one or more of the following US patents:
6,510,007; 7,787,197; 6,683,314; 6,809,804; 7,129,505; 6,897,954; 8,528,427


 Becton, Dickinson and Company, BD Biosciences, 2350 Qume Drive, San Jose, CA 95131 USA

 Becton Dickinson Ireland Ltd., Donore Road, Drogheda, Co. Louth, A92 YW26, Ireland

 BD Switzerland Sàrl, Terre Bonne Park – A4, Route de Crassier 17, 1262 Eysins, Switzerland

BD, the BD Logo, and FACSLytic are trademarks of Becton, Dickinson and Company or its affiliates.
© 2021 BD. All rights reserved.

Made in USA
23-20838-02



(01)00382906511652(11)221123(21)R651165000266

Макет заводской маркировки Цитофлуориметра проточного BD FACSLytic
(пример для варианта исполнения: «Цитофлуориметр проточный BD FACSLytic
(2 лазера, 6 флуоресцентных параметра)»)

Model:  **BD FACSLytic™ Flow Cytometer**

REF 664587

Complies with FDA performance standards for laser products except for deviations pursuant to laser notice No. 50, dated June 24, 2007. Class 1 Laser Product per IEC/EN 60825-1:2014. This device complies with Part 15 of the FCC Rules. Operation is subjected to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This Class A digital apparatus meets all requirements of the Canadian Interference Causing Equipment Regulations.

Appareil à laser de classe 1, IEC/EN 60825-1:2014. Cet appareil numérique de la classe A respecte toutes les exigences du Règlement sur le matériel brouilleur du Canada.



SN:

2022-12-06

R654587000243

VAC-Hz:

100-240 VAC (50-60 Hz)

Total Power:

200 W



The BD FACSLytic™ flow cytometer is covered by one or more of the following US patents:
6,510,007; 7,787,197; 6,683,314; 6,809,804; 7,129,505; 6,897,954; 8,528,427


 **Becton, Dickinson and Company, BD Biosciences, 2350 Qume Drive, San Jose, CA 95131 USA**

 **Becton Dickinson Ireland Ltd., Donore Road, Drogheda, Co. Louth, A92 YW26, Ireland**

 **BD Switzerland Sàrl, Terre Bonne Park – A4, Route de Crassier 17, 1262 Eysins, Switzerland**

BD, the BD Logo, and FACSLytic are trademarks of Becton, Dickinson and Company or its affiliates.
© 2021 BD. All rights reserved.

Made in USA
23-20838-02



(01)00382906545879(11)221206(21)R654587000243

Макет заводской маркировки Цитофлуориметра проточного BD FACSLytic
(пример для варианта исполнения: «Цитофлуориметр проточный BD FACSLytic
(3 лазера, 8 флуоресцентных параметра)»)

Model:  **BD FACSLytic™ Flow Cytometer****REF 659180**

Complies with FDA performance standards for laser products except for deviations pursuant to laser notice No. 50, dated June 24, 2007. Class 1 Laser Product per IEC/EN 60825-1:2014. This device complies with Part 15 of the FCC Rules. Operation is subjected to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This Class A digital apparatus meets all requirements of the Canadian Interference Causing Equipment Regulations.

Appareil à laser de classe 1, IEC/EN 60825-1:2014. Cet appareil numérique de la classe A respecte toutes les exigences du Règlement sur le matériel brouilleur du Canada.

**2022-11-29**

SN:

R659180000762

VAC-Hz:

100–240 VAC (50–60 Hz)

Total Power:

200 W



The BD FACSLytic™ flow cytometer is covered by one or more of the following US patents:
6,510,007; 7,787,197; 6,683,314; 6,809,804; 7,129,505; 6,897,954; 8,528,427



Becton, Dickinson and Company, BD Biosciences, 2350 Quime Drive, San Jose, CA 95131 USA




Becton Dickinson Ireland Ltd., Donore Road, Drogheda, Co. Louth, A92 YW26, Ireland



BD Switzerland Sàrl, Terrie Bonine Park – A4, Route de Crassier 17, 1262 Eysins, Switzerland

BD, the BD Logo, and FACSLytic are trademarks of Becton, Dickinson and Company or its affiliates.
© 2021 BD. All rights reserved.

Made in USA
23-20838-02

**(01)00382906591807(11)221129(21)R659180000762**

Макет заводской маркировки Цитофлуориметра проточного BD FACSLytic
(пример для варианта исполнения: «Цитофлуориметр проточный BD FACSLytic
(3 лазера, 10 флуоресцентных параметра)»)

Model:  **BD FACSLytic™ Flow Cytometer**

REF 663029

Complies with FDA performance standards for laser products except for deviations pursuant to laser notice No. 50, dated June 24, 2007. Class 1 Laser Product per IEC/EN 60825-1:2014. This device complies with Part 15 of the FCC Rules. Operation is subjected to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This Class A digital apparatus meets all requirements of the Canadian Interference Causing Equipment Regulations.

Appareil à laser de classe 1, IEC/EN 60825-1:2014. Cet appareil numérique de la classe A respecte toutes les exigences du Règlement sur le matériel brouilleur du Canada.



2022-11-30

SN:

R663029000666

VAC-Hz:

100-240 VAC (50-60 Hz)

Total Power:

200 W



eifu.bd.com



The BD FACSLytic™ flow cytometer is covered by one or more of the following US patents: 6,510,007; 7,787,197; 6,683,314; 6,809,804; 7,129,505; 6,897,954; 8,528,427


 **Becton, Dickinson and Company, BD Biosciences, 2350 Qume Drive, San Jose, CA 95131 USA**

 **Becton Dickinson Ireland Ltd., Donore Road, Drogheda, Co. Louth, A92 YW26, Ireland**

 **BD Switzerland Sàrl, Terre Bonne Park - A4, Route de Crassier 17, 1262 Eysins, Switzerland**

BD, the BD Logo, and FACSLytic are trademarks of Becton, Dickinson and Company or its affiliates.
© 2021 BD. All rights reserved.

Made in USA
23-20838-02



(01)00382906630292(11)221130(21)R663029000666

Макет заводской маркировки Цитофлуориметра проточного BD FACSLytic
(пример для варианта исполнения: «Цитофлуориметр проточный BD FACSLytic
(3 лазера, 12 флуоресцентных параметра)»)

Model:  **BD FACS™ Universal Loader**

REF 651166

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This Class A digital apparatus meets all requirements of the Canadian Interference Causing Equipment Regulations.

Cet appareil numérique de la classe A respecte toutes les exigences du règlement sur le matériel brouilleur du Canada.



SN:

2022-12-13


R651166001693

Voltage:

12 VDC
24 VDC

Power:

18 W
36 W



Becton, Dickinson and Company, BD Biosciences, 2350 Qume Drive, San Jose, CA 95131 USA



Becton Dickinson Ireland Ltd., Danvers Road, Drogheda, Co. Louth, A92 YW26, Ireland



BD Switzerland Sàrl, Terre Barne Park - AA, Route de Crassier 17, 1262 Eysins, Switzerland

BD, the BD Logo and FACS are trademarks of Becton, Dickinson and Company or its affiliates.
© 2021 BD. All rights reserved.

Made in USA
23-23986-01



(01)00382906511669(11)221213(21)R651166001693

Макет заводской маркировки модуля для подачи проб универсального BD FACS Universal Loader (Пример)

Модуль для подачи проб универсальный BD FACS Universal Loader

Производитель:

Becton, Dickinson and Company, BD Biosciences
(Бектон, Дикинсон энд Компани, БД Байосайенсиз)
2350 Qume Drive, San Jose, California, 95131, США

Страна происхождения:

Сделано в США

Номер по каталогу **REF**: см. на заводской маркировке

Серийный номер **SN**: см. на заводской маркировке

Условия эксплуатации:

Температура: от 15°C до 30°C;
Относительная влажность:
от 15% до 85% (без конденсации);
Атмосферное давление от 0,8 до 1,0 амт.

Условия хранения без эксплуатации:

Температура: от 0°C до 49°C;
Относительная влажность:
от 10% до 20 % (без конденсации).



Уполномоченный представитель производителя:

ООО «Бектон Дикинсон Восток»
127051 Российская Федерация, Москва, ул. Садовая-Самотечная, 24/27,
2 этаж, помещение I, комната 18
Тел.: +7 (495)775-85-82; e-mail: complaints_CIS@bd.com

Регистрационное Удостоверение № _____ от ДД.ММ.ГГГГ

Информацию на русскоязычном стикере просим считать приоритетной при поставке в РФ

XXXXXX

xxxx&xxxx

Макет русскоязычной маркировки: компонента медицинского изделия
«Модуль для подачи проб универсальный BD FACS Universal Loader»

ПРОШТО, ПРОЧУМІРОВАНО И
СКРЕПЛЕНО ПОДПИСЬЮ И ПЕЧАТЬЮ
В. Листов
СПЕЦИАЛИСТ ПО РЕГИСТРАЦИИ МИ
Вилкова Д.И.
МОСКВА * 502470091101