

УТВЕРЖДАЮ

Генеральный директор
ООО «Компания «БиВи»

Т.А. Носова

«24» апреля 2023 г.



по ДОВЕРЕННОСТИ
77 АД 1754197

КУСАЙЛО Р. В.

Фотографические изображения изделия:

**Реагент лизирующий для анализатора автоматического гематологического 5-Diff для
диагностики in vitro (Hematology analyzer lyse (LYA-3/ LYR-3))**

Shenzhen Dymind Biotechnology Co., Ltd., («Шэньчжэнь Дайминд Биотекнологджи Ко., Лтд.»),
Китай



DYMIND

500 mL

LOT 2022052401

20220524



20240523

LYA-3



WARNINGS:

1. Avoid skin or eye contact. In case of skin or eye contact, flush the affected area with plenty of water immediately and go see a doctor. In case of inhaling, immediately go see a doctor.
2. Finish product before the expiration date. Once it is opened, it should be finished within 60 days.
3. If product is frozen (either partially or completely), thaw completely, warm to room temperature and mix thoroughly by gentle inversion. Verify background results before analyzing patient samples.
4. Dispose of waste product in accordance with governmental requirements

EC REP

Eunitor GmbH

Kennedydamm 5, 40476 Duesseldorf, Germany



IVD



+2°C



+30°C

MADE IN CHINA

Реагент лизирующий LYA-3, 500 мл

Hematology Analyzer Lyse

LYA-3 Lyse is applicable to the DH5X series, DH7X series, HS-5, HS-5CRP, HS-7, DX-CRP series, DHXXCRP series, UN7X series hematology analyzers manufactured by Dymind. Among which, X represents the digit and/or letter. It is used to break down red blood cell and cooperate with the LYA-2 Lyse to 4-differentiate WBCs. Use the LYA-3 Lyse as instructed by the operation manual of the analyzer.

ACTIVE INGREDIENTS:

Surface active agent; Buffering agent; Antifungal and antibacterial Agents.

STORAGE AND STABILITY:

The product will be stable for 2 years when stored at 2°C to 30°C (35°F to 86°F).

If the operating environment is 15°C to 30°C (59°F to 86°F), discard opened container after 60 days.




Shenzhen Dymind Biotechnology Co., Ltd.

10th Floor, Building B, High-tech Park, Guangqiao Road, Tianliao Community
Yutang Street, Guangming District, Shenzhen 518107, P. R. China

Реагент лизирующий LYA-3, 500 мл



Реагент лизирующий LYA-3, 500 мл



Реагент лизирующий LYA-3, 500 мл



DYMIND

1 L

LOT 2022112301



20221123



20241122

LYA-3




WARNINGS:

1. Avoid skin or eye contact. In case of skin or eye contact, flush the affected area with plenty of water immediately and go see a doctor. In case of inhaling, immediately go see a doctor.
2. Finish product before the expiration date. Once it is opened, it should be finished within 60 days.
3. If product is frozen (either partially or completely), thaw completely, warm to room temperature and mix thoroughly by gentle inversion. Verify background results before analyzing patient samples.
4. Dispose of waste product in accordance with governmental requirements

EC REP

Eunitor GmbH

Kennedydamm 5, 40476 Duesseldorf, Germany



MADE IN CHINA

Реагент лизирующий LYA-3, 1 л

Hematology Analyzer Lyse

LYA-3 Lyse is applicable to the DHSX series, DH7X series, HS-5, HS-5CRP, HS-7, DX-CRP series, DHXXCRP series, UN7X series hematology analyzers manufactured by Dymind. Among which, X represents the digit and/or letter. It is used to break down red blood cell and cooperate with the LYA-2 Lyse to 4-differentiate WBCs. Use the LYA-3 Lyse as instructed by the operation manual of the analyzer.

ACTIVE INGREDIENTS:

Surface active agent; Buffering agent; Antifungal and antibacterial Agents.

STORAGE AND STABILITY:

The product will be stable for 2 years when stored at 2°C to 30°C (35°F to 86°F).


If the operating environment is 15°C to 30°C (59°F to 86°F), discard opened container after 60 days.



Shenzhen Dymind Biotechnology Co., Ltd.

10th Floor, Building B, High-tech Park, Guangqiao Road, Tianliao Community
Yutang Street, Guangming District, Shenzhen 518107, P. R. China

Реагент лизирующий LYA-3, 1 л



Реагент лизирующий LYA-3, 1 л



Реагент лизирующий LYA-3, 1 л



DYMIND

500 mL

LOT

2022052402



20240523

LYR-3



20220524




WARNINGS:

1. Avoid skin or eye contact. In case of skin or eye contact, flush the affected area with plenty of water immediately and go see a doctor. In case of inhaling, immediately go see a doctor.
2. Finish product before the expiration date. Once it is opened, it should be finished within 60 days.
3. If product is frozen (either partially or completely), thaw completely, warm to room temperature and mix thoroughly by gentle inversion. Verify background results before analyzing patient samples.
4. Dispose of waste product in accordance with governmental requirements.

EC REP

Eunitor GmbH

Kennedydamm 5, 40476 Duesseldorf, Germany



MADE IN CHINA

Реагент лизирующий LYR-3, 500 мл

Hematology Analyzer Lyse

LYR-3 Lyse is applicable to the DM79, DM78, DM72, DM71, DM79X, DM78X, DM72X, DM71X hematology analyzers manufactured by Dymind. It is used to break down red blood cell and cooperate with the LYR-2 Lyse to 4-differentiate WBCs. Use the LYR-3 Lyse as instructed by the operation manual of the analyzer.

ACTIVE INGREDIENTS:

Quaternary ammonium surfactants and NaCl.

STORAGE AND STABILITY:

The product will be stable for 2 years when stored at 2°C to 30°C (35°F to 86°F).

If the operating environment is 15°C to 30°C (59°F to 86°F), discard opened container after 60 days.



Shenzhen Dymind Biotechnology Co., Ltd.

10th Floor, Building B, High-tech Park, Guangqiao Road, Tianliao Community
Yutang Street, Guangming District, Shenzhen 518107, P. R. China

Реагент лизирующий LYR-3, 500 мл



Реагент лизирующий LYR-3, 500 мл



Реагент лизирующий LYR-3, 500 мл



Реагент лизирующий LYR-3, 1 л

Hematology Analyzer Lyse

LYR-3 Lyse is applicable to the DM79, DM78, DM72, DM71, DM79X, DM78X, DM72X, DM71X hematology analyzers manufactured by Dymind. It is used to break down red blood cell and cooperate with the LYR-2 Lyse to 4-differentiate WBCs. Use the LYR-3 Lyse as instructed by the operation manual of the analyzer.

ACTIVE INGREDIENTS:

Quaternary ammonium surfactants and NaCl.

STORAGE AND STABILITY:

The product will be stable for 2 years when stored at 2°C to 30°C (35°F to 86°F).

If the operating environment is 15°C to 30°C (59°F to 86°F), discard opened container after 60 days.



Shenzhen Dymind Biotechnology Co., Ltd.

10th Floor, Building B, High-tech Park, Guangqiao Road, Tianliao Community,
Yutang Street, Guangming District, Shenzhen 518107, P. R. China

Реагент лизирующий LYR-3, 1 л




Реагент лизирующий LYR-3, 1 л



Реагент лизирующий LYR-3, 1 л


Hematology Analyzer Lyse



List	Model	Specification	LOT	Production date	Expiry date	PCS
	LYR-3	500mL	2022052401	24.05.2022	23.05.2024	3
	LYR-3	500mL	2022052402	24.05.2022	23.05.2024	7



Made in China



Eunitor GmbH
Kennedydamm 5, 40476 Duesseldorf, Germany



Shenzhen Dymind Biotechnology Co., Ltd.


10th Floor, Building B, High-tech Park, Guangqiao Road, Tianliao Community, Yutang Street,
Guangming District, Shenzhen 518107, P.R.China

Country : Russia-Company BW, LLC



No. : DM220606023-1

Транспортная маркировка «Реагент лизирующий для анализатора автоматического гематологического 5-Diff для диагностики in vitro (Hematology analyzer lyse (LYA-3/ LYR-3))», вариант исполнения Реагент лизирующий LYR-3, 500 мл

Hematology Analyzer Lyse



List	Model	Specification	LOT	Production date	Expiry date	PCS
	LYA-3	500mL	2022052401	24.05.2022	23.05.2024	3
	LYA-3	500mL	2022052402	24.05.2022	23.05.2024	7



Made in China



Eunitor GmbH
Kennedydamm 5, 40476 Duesseldorf, Germany



Shenzhen Dymind Biotechnology Co., Ltd.

10th Floor, Building B, High-tech Park, Guangqiao Road, Tianliao Community, Yutang Street,
Guangming District, Shenzhen 518107, P.R.China

Country : Russia-Company BW, LLC

No. : DM220606023-1

C/N : 78-71

Транспортная маркировка «Реагент лизирующий для анализатора автоматического гематологического 5-Diff для диагностики in vitro (Hematology analyzer lyse (LYA-3/ LYR-3))», вариант исполнения Реагент лизирующий LYA-3, 500 мл

Всего прошито, пронумеровано и скреплено
печатью и подписью: 19 листа (-ов)

Генеральный директор _____ / Носова Т.А
М.П.

ПО ДОВЕРЕННОСТИ
77 АД 1754197
КУСНИКОВ В.

