







система
фиксации
(часть индивидуального
уровня)


стенка

phenox

Manufactured by
phenox GmbH
Lise-Meitner-Allee 31
44801 Bochum, Germany
Fon: +49 (0) 234/36919-0
Fax: +49 (0) 234/36919-19
www.phenox.info

p64 Flow Modulation Device

P64-450-27



Delivery System with distal wire and tip

Use in

Vessel Ø (mm)	Length p64 (mm)
4.50 (max. Vessel Ø)	27
4.25	33
4.00 (min. Vessel Ø)	40

Other Conditions

p64 Ø (mm)	Length p64 (mm)
4.8 (Unconstrained)	13.2
0.69 (in ID 0.027 inch Microcatheter)	71.4

BB10 D - 2013/07

REF	LOT	CONT		MICROCATHETER COMPATIBILITY
P64-450-27	307349	1	2016-07	ID 0.027"

STERILE | EO

NONPYROGENIC

CE 0297



22°C



(01)04260123780470(17)160731(10)307349

phenox

REF P64-450-27



p64

LOT

307349



2016-07

сгуст

узкая часть
губчатая

изгибывание графика

phenox

Manufactured by
phenox GmbH
Lindenstrasse 11
44891 Bochum, Germany
Fax: +49 (0) 234 38919-0
Fax: +49 (0) 234 38919-10
www.phenox.info

p64 Flow Modulation Device

P64-450-27

Delivery System with distal wire and tip

Use in	
Vessel Ø (mm)	Length p64 (mm)
4.5	27
5.5	33
6.5	40

Other Conditions	
p64 Ø (mm)	Length p64 (mm)
4.5	13.2
5.5	13.2
6.5	13.2

REF	LDT	CONF	MONOCATHETER COMPATIBILITY
P64-450-27	307349	1	2016-07

STERILE EO 0297 NONPYROGENIC

22°C

phenox P64-450-27 p64 307349 2016-07

p64 Flow Modulation Device p64 P64-450-27 307349

p64 Flow Modulation Device p64 P64-450-27 307349


p64 Flow Modulation Device p64 P64-450-27 307349

p64 Flow Modulation Device p64 P64-450-27 307349

p64 Flow Modulation Device p64 P64-450-27 307349


p64 Flow Modulation Device p64 P64-450-27 307349

p64 Flow Modulation Device p64 P64-450-27 307349



сгуст

состав
продукта
изготовленного




phenox

Manufactured by
phenox GmbH
Lisa-Meiner-Allee 31
44801 Bochum, Germany
Phone: +49 (0) 234/36919-0
Fax: +49 (0) 234/36919-19
www.phenox.info

p64 Flow Modulation Device

P64-450-27



Delivery System with distal wire and tip

Use in	
Vessel Ø (mm)	Length p64 (mm)
4.50 (max. Vessel Ø)	27
4.25	33
4.00 (min. Vessel Ø)	40

Other Conditions

p64 Ø (mm)	Length p64 (mm)
4.8 (Unconstrained)	13.2
0.69 (in ID 0.027 inch Microcatheter)	71.4

REF	LOT	CONT	DATE	MICROCATHETER COMPATIBILITY
P64-450-27	307349	1	2016-07	ID 0.037"

STERILE | EO 0297 22°C
NONPYROGENIC



(01)04260123780470(17)160731(10)307349

phenox

REF P64-450-27

p64

p64 Flow Modulation Device p64



(01)042601237

.50-27 LOT 307349

p64 Flow Modulation Device p64



(01)042601237

P64-450-27 LOT 307349

p64 Flow Modulation Device p64



(01)042601237

P64-450-27 LOT 307349


p64 Flow Modulation Device p64



(01)042601237

P64-450-27 LOT 307349

p64 Flow Modulation Device p64



(01)042601237

P64-450-27 LOT 307349

p64 Flow Modulation Device p64



(01)042601237

P64-450-27 LOT 307349

Результаты, полученные в
испытании верности: результаты
в норме.
Инженер *[подпись]*
ООО, дата *[подпись]*

